

Validation of Simple Sequence Length Polymorphism (SSLP) regions of commonly used mouse strains for marker assisted speed congenics screening

Channabasavaiah B Gurumurthy^{1,4*}, Poonam S Joshi^{1,4}, Scott Kurz¹, Masato Ohtsuka², Rolen M Quadros¹, Donald W Harms¹, Kent Lloyd³

(1) Mouse Genome Engineering Core Facility, Department of Genetics, Cell Biology and Anatomy, University of Nebraska Medical Center, Omaha, NE, USA

(2) Department of Molecular Life Science, Division of Basic Medical Science and Molecular Medicine, Tokai University School of Medicine, 143 Shimokasuya, Isehara, Kanagawa 259-1193, Japan

(3) Mouse Biology Program, University of California, Davis, 2795 2nd Street, Davis, CA, 95618, USA

(4) These authors contributed equally

* Corresponding author

Running Title: SSLP markers panel for common mouse strains

Key Words: SSLP, Speed Congenic, Microsatellite Marker, Mouse Strain

Abstract

Marker-assisted speed congenics technique is commonly used to facilitate backcrossing of mouse strains in nearly half the time it normally takes otherwise. Traditionally, the technique is performed by analyzing PCR amplified regions of Simple Sequence Length Polymorphism (SSLP) markers between the recipient and donor strains: **offspring with the highest number of markers showing the recipient genome across all chromosomes is chosen for the next generation.**

Although there are well defined panels of SSLP makers established between certain pairs of mice strains, they are incomplete for most strains. The availability of well- established marker sets for speed congenic screens would enable the scientific community to transfer mutations across strain backgrounds. In this study, we tested the suitability of over 400 SSLP marker sets among 10 mouse strains commonly used for generating genetically engineered models. The panel of markers presented here can readily identify the specified strains and will be quite useful in marker assisted speed congenic screens. Moreover, unlike newer Single Nucleotide Polymorphism (SNP) array methods which require sophisticated equipment, the SSLP markers panel described here only uses PCR and agarose gel electrophoresis of amplified products; therefore it can be performed in most research laboratories.

Introduction

In the recent years, there has been a steady increase in the creation and use of genetically engineered mutant mice for use in biomedical research. Frequently, the genetic background of such mice allow only specific experiments and the mutations have to be transferred to different genetic background(s) to facilitate other kinds of experiments. There are many examples where genetic background is shown to influence the phenotype of a transgenic or knockout mouse line (Cao et al., 2013; Davie et al., 2007; Doetschman, 2009; Fleming et al., 2001; Gale et al., 2009; Heiman-Patterson et al., 2011; Johnson et al., 2006; MacLennan et al., 2011; Montagutelli, 2000; Puccini et al., 2013; Siuciak et al., 2008; Thifault et al., 2002). Traditionally however, mutant mice are generated using strains that have shown exceptional performance in terms of their suitability for production of transgenic or knockout mice lines. For example, FVB strain and BDF1 strain mice are most commonly used for transgenic mice production (Auerbach et al., 2003) and ES cells derived from 129 inbred strains are commonly used for knockout mice production. For technical reasons, chimeras developed in knockout mice generation will carry a mixed genetic background (e.g. 129 and B6) adding further complexity to the analysis (Glaser et al., 2005). Furthermore, ES cells derived from C57BL6/N inbred strain have been used in mouse genetic resources such as KOMP (Knock-Out Mouse Project) and EUCOMM (European Conditional Mouse Mutagenesis)(Beier, 2010). Also, particular strains of mice seem to have better suitability for a specific research purpose. In addition, mixed background strains between C57BL/6J and DBA/2J were used for ENU mutagenesis projects because the cross between two different inbred strain (one strain [male] was used for ENU exposure and mated with another strain [female]) is useful for mapping and positional cloning of the mutated gene using genetic polymorphisms existing between them (Inoue et al., 2004; Masuya et al., 2005) The main disadvantage of a given mutation under a particular strain background is that the strain background may limit its use for a specific research purpose. In such a situation the mutation is needed to be transferred into a strain background of choice through a process called backcrossing.

Backcrossing involves about 10 generations of successive breeding into a recipient strain of choice to achieve 99.9% congenic (genetic composition) for that strain (www.informatics.jax.org/silverbook/). This painstaking process consumes about 2.5 to 3 years of time; a fact that often limits its feasibility and usefulness given the pace of scientific research. In some cases, studies are published with animals after only 5 generations of backcrossing, in an attempt to compensate the time required and the need to obtain some results in the new strain (Cendán et al., 2005). A technique called “Marker assisted speed congenic” used since over a decade helps achieving congenic strain in 5, or less, unlike the usual 10 generations required in traditional backcrossing (Markel et al., 1997; Wakeland et al., 1997). The small sequence differences between mouse strains called ‘microsatellite markers’ served as useful tools in detecting the chromosome regions of origin in the offspring when two inbred strains of mice are bred together. To use in these assays, many microsatellite markers have been identified and characterized by various researchers between the donor and recipient strains of their choice (Collins et al., 2003; Estill and Garcia, 2000, p. 6; Farber et al., 2006; Jakoubova et al., 2000; Ogonuki et al., 2009; Sakai et al., 2004; Suemizu et al., 2008; Teppner et al., 2004; Witmer et al., 2003; Lamacchia et al., 2007). However, the information about the markers that can differentiate between commonly used strains for transgenic and knockout mice generation is not tested sufficiently and is not available readily. In this study, we tested 423 markers, ~ 10 to 30 per chromosome using the genomic DNAs from 10 commonly used mice strains—particularly the strains used in transgenic and knock-out research. We evaluated the markers that could be used in agarose gel electrophoresis method which is a simple technique commonly used in most molecular biology laboratories these days. The data presented here will serve as a valuable tool for various investigators in choosing the markers useful for their speed congenic breeding.

Materials and methods

Selection of Oligonucleotide primers: The primers were chosen based on the following criteria. 1) To evaluate and establish at least 6 markers per chromosome, 2) distance between the adjacent markers was kept as minimum as 10 to 15 centi Morgans (cM). 3) Polymorphic bands appreciable when analyzed in a 4% agarose gel electrophoresis and resolved by electrophoresis distance of up to 10 centimeters from the loading wells. All markers were chosen from the Mouse Genome Informatics (MGI) database links [http://www.informatics.jax.org/searches/probe_report.cgi?_Refs_key=22816 and <ftp://ftp.informatics.jax.org/pub/datasets/index.html> (numbers 7 and 12 in the list)].

Mice strains: The mice were purchased from Charles River or The Jackson Laboratories. The strains, the rationale for including these strains in the panel and the vendors are listed in the table below.

Strain	Vendor, Stock Number	Rationale
C57BL/6J	JaxMice, 000664	Commonly used strain in biomedical research
FVB/NJ	JaxMice, 001800	Commonly used strain for transgenic mice generation
129 X1/SvJ	JaxMice, 000691	Strains from which highly successful germ line transmission competent ES cells were derived.
129S2	Charles River, 476	
DBA2/J	JaxMice, 000671	Together with C57BL6, these strains are used as hybrid strains for transgenic mouse production
CBA/J	JaxMice, 000656	
SJL/J	JaxMice, 000686	
Balb/c	Charles River, 028	Is among the top 2-3 most widely used inbred strain
C3H/HeJ	Charles River 025	Used in a variety of research areas including cancer, infectious disease and cardiovascular biology.
NOD	Charles River	

DNA extraction, PCR amplification and agarose gel electrophoresis

The DNA samples were extracted from the tail pieces of about 3-5 mm long using Gentra Puregene Tissue Kit (Cat. # 158622). PCR was performed in a volume of 20 μ l 1X reaction buffer containing 20 mM Tris pH 8.4, 50 mM KCl, 3 mM MgCl₂ and 1 unit of Taq DNA polymerase (New England Biolabs, Catalogue Number M0273) with the following conditions, one cycle of 95°C –2 min followed by 35 cycles of 95°C –30 sec, 55°C – 30 sec, 72°C–60 sec and one cycle of 72°C – 5 minutes, followed by a holding temperature at 4°C until the samples were removed from the machine. Fifteen μ l of PCR products were resolved using 4% agarose gels for 120 to 150 minutes at 200 constant volts electric current. The agarose was purchased from Phenix Research Products (Item Number. RBA-500: Molecular Biology grade) and the gels were prepared on 0.5X TAE buffer diluted from a stock of 50X (Fisher Scientific, Cat. Numner. BP1332-20). The gels contained Ethidium Bromide dye (0.5 μ g/ml) to aid in visualization of PCR bands. Each panel of gel included one or more lanes of 100 base pair Molecular Weight marker (New England Biolabs, Catalogue Number N3231) to assess the PCR product sizes. The bromophenol blue dye-front was let to run for up to 10 centimeters from the wells and the gels were imaged using BioRad Gel Doc XR system. Wherever necessary, the gels were run longer to resolve the bands.

Analysis and interpretation of polymorphic PCR bands

The cropped images were imported into an excel file for analysis and interpretation. Band sizes were assigned numbers 1, 2, 3 or 4 to indicate their sizes compared to the rest of the bands in that set. Number 1 was assigned to the smallest sized polymorphic band, 2 to the next biggest in the group and so on. The excel file along with the gel images was converted into a .pdf file for generating the Figure 1.

Results and discussion

Mutant mice created using transgenic and knockout techniques are available under certain specific backgrounds. In order to best use such mutants for a specific research purpose, they routinely need to be bred into other strain backgrounds through successive breeding of about 10 generations. A quicker way to attain highest recipient genome can be achieved by a process called speed congenic breeding in which the polymorphic markers between the recipient and donor strains are screened among offspring in each generation and the offspring with highest recipient genome is chosen as a breeder for the subsequent generation (Collins et al., 2003; Markel et al., 1997; Wakeland et al., 1997). Although there are a few reports describing the marker sets suitable for a certain pairs of strains, there are no well-established marker sets available across the most commonly used strains in transgenic and knockout mouse techniques. Our primary objective of this work was to test the suitability of several markers in an agarose gel electrophoresis method; a technique that uses least expensive equipment, reagents and is readily available in most molecular biology laboratories. We sought to test a large number of microsatellite markers (Tables 1, 3 and S1) among 10 most commonly used mouse strains particularly those used in transgenic and knockout mice techniques.

We used Mouse Genome Informatics (MGI) database and short listed the primers based on the criteria described in Materials and Methods. MGI database lists several of the microsatellite markers identified in various inbred strains. Although the database is extensive, it is difficult to choose marker sets for speed congenic screening because it lacks the critical information such as the sizes of PCR products of various markers in different inbred strains and if the differences among strains can be identified using conventional agarose gel electrophoresis.

As per the information available on MGI database, predicted PCR band size differences among strains for some markers ranged from as less as a few base pairs to over 100 base pairs. Selection criteria for markers and the mice strain is described in Materials and Methods section. Taking C57BL6J and 129X1 strains as a comparison pair for example, we chose markers with differences of 30 bp and above (as per MGI database), a range that can be easily resolved using about 2% agarose gels. We aimed to choose markers in this range for all the chromosomal locations with an interval of 5 to 15 cM. However, we were unable to find suitable markers in some regions with this criterion. In such locations, we tested markers with as less as 8 to 12 bp size difference. Such small differences can be best resolved using poly acrylamide gel electrophoresis (PAGE). However, in a typical speed congenic screening that involves analysis of about 100 markers for each sample, applying PAGE method to screen becomes quite laborious. We used 4% agarose gels to resolve markers even with very small size differences; this seemed to be sufficient to resolve the bands when they were run about 8 to 10 centimeters from the loading well. In order to keep the assay conditions uniform we used 4% agarose gel for all the markers tested.

Table 1 shows the list of the markers tested and found polymorphic in at-least one of the 10 strains analyzed. We tested a total of 423 markers of which we present useful data for 195 markers. The gel images and interpretation of polymorphisms are tabulated in Figure 1. The gel images indicate that there were readily appreciable size differences for most markers whereas some markers showed very narrow size difference between some strains. In general, sizes of the majority of the markers matched the information on MGI database but there were some discrepancies noticed. The agarose gel electrophoresis using 4% gels could be useful in detecting differences among the markers in different strains.

The data presented in Figure 1 indicate that some markers could readily distinguish several strains from each other. We initially aimed to identify good panel of markers for every 10 cM. Although we found adjacent markers as close as 10 cM in most locations, we failed to achieve this in some regions. We screened several available markers in such regions and were unable to find markers that matched the criteria set forth in our assay. It should be noted that there were markers that differed very minutely among some strains. If such markers are chosen in an actual speed congenic screening, we recommend including samples from both donor and recipient strains and an equimolar mixture of the two in the panel in order to ensure the proper banding pattern of offspring. The gel images and the interpretation of polymorphism provided in figure 1 will serve as a comprehensive tool for a researcher who wishes to undertake congenic breeding between a pair of strains from the list.

The data from the markers that did not meet our criteria fell into one of the following categories: 1) Yielded same size bands in all the strains. 2) Failed to amplify a product in some strains. 3) Amplified multiple PCR products. 4) **Did not yield a reliable PCR product in most/all of the strains.** The category 1 means that the bands were not resolvable by 4% agarose gels. We presume that such markers would be < 6 base-pairs different from one another (Neuhaus et al 1997 Mammalian Genome 8, 506-509 (1997)) or they do not differ in size at all. Category 2 would probably mean that the primers did not have perfect binding sites in those strains where there was no PCR product. Category 3 makers may need further optimization of PCR conditions. However, one of our aims was to keep the conditions uniform to all the markers in order to keep the assay simple: it makes the assay cumbersome if PCR conditions have to be varied among different marker sets. It should be noted that the category 4 markers (that did not **yield a reliable PCR product in most/all of the strains**) **might have not worked due to technical errors that** may have occurred during some steps like primer synthesis or primer reconstitution; we do not rule out the possibility that some of these primers may yield PCR products if they are re-synthesized. The markers that did not meet our criteria because of the above reasons are listed in **Table 2 and** Supplementary Table S1: We felt that providing the list of primers and markers that failed in our hands will help researchers to compare the results if they intend to test more markers to expand the streamlined panel presented here.

There is some debate about the minimum number of markers needed per chromosome for identifying the regions between the donor and recipient strains (Armstrong et al., 2006). A study concluded that, as less as three markers per chromosome were sufficient to achieve similarly meaningful results as that of over 6 markers per chromosome (Goto et al., 2005). Computer simulation (Wakeland et al., 1997) indicated that a relatively modest selection effort of 60 evenly spaced markers with 25cM spacing (corresponds to 3 markers per chromosome), 16 males per generation would typically reduce unlinked donor genome contamination to below 1% by four backcross generations (N5). We conclude that the list presented here can serve to choose panels of markers for most two strain combinations with at least 3 to 4 markers per chromosome (with exception of a very few combinations). Further studies that compare smaller and larger panel of markers between the same set of samples for marker assisted speed congenics are needed to address this question unequivocally.

In the recent years there have been advancements in the approaches used for speed congenics. These methods include i) use of fluorescently labeled primers to amplify SSLP markers followed by resolving the products in sequencing gels (Witmer et al., 2003) and ii) microarray chips of SNPs (single Nucleotide Polymorphisms) (Petkov, 2004; Petkov et al., 2004; Tsang et al., 2005; Witmer et al., 2003). With the advent of newer methods particularly those that use SNP based marker analysis, very high number of markers per chromosome can be screened simultaneously that increases the resolution by several fold compared to the conventional SSLP based markers analysis. Although there are some computer simulation studies for assessing the efficiency of speed congenic screening in general (Armstrong et al.,

2006) and algorithm based reports to compare SNP and SSLP based speed congenic screens(Gorham et al., 2012) there are no systematic studies to compare the two methods to assess the efficiency and cost effectiveness of each method. Here, we compare the SNP and SSLP based approaches in terms of their adoptability and feasibility to most laboratory settings including their cost. The hands-on-time in performing gel based assays has been reduced greatly by newer methods that use SNP arrays. However these methods have some limitations compared to traditional agarose gel electrophoresis (AGE) based SSLP marker analysis. 1) The newer methods are expensive in terms of the initial investment in reagents costs and/or operational costs compared to SSLP-AGE method. On the contrary, basic requirements needed for AGE based systems are readily available in any molecular biology laboratory and the only additional investment needed will be to synthesize the required oligonucleotides 2) subset of markers to be analyzed in the subsequent generations cannot be skipped from the panel unlike in the AGE based method where the number of markers to be analyzed will become significantly reduced in successive generations and so also the overall cost of the assay. Other advantages of AGE based SSLP marker analysis over these newer methods are i) the equipment needed to perform the assay is readily available in most laboratories and ii) can be routinely performed by many researchers and technicians without the need of special training as needed for SNP based approaches.

Considering the highest resolution that is possible with the SNP based method, it can be regarded as the superior method of all. However, the microchips that are currently available are expensive and the cost for analyzing each mouse DNA sample runs about US \$100 to \$150. Assuming about analysis of 15 mouse DNA samples per generation for 5 generations, a typical speed congenic project would need about \$15,000 to \$22,500. On the other hand, when using SSLP-AGE approach, since marker analysis is done manually, the markers that were fixed in the previous generation can be skipped in the successive generations; this makes the SSLP-AGE based method cost effective compared to SNP based method. It is estimated that, a typical SSLP based speed congenic screen that employs analysis of 15 samples per generation for 5 generations using 5 markers per chromosome would need about 2000 to 2300 PCR reactions. At the rate of \$2 to \$2.5 per reaction (cost analysis done in our laboratory) will cost about \$4000 to \$5750 for one full speed congenic project. Furthermore, SSLP-AGE method can be performed in any simple molecular biology labs compared to SNP method that requires expensive equipments.

Conclusions: Although, some information about microsatellite marker differences between the commonly used inbred mouse strains were available, there was no systematic study to validate a large panel of markers for SSLP-AGE based speed congenic screening. The panel of marker sets validated and presented in this study serves as a ready reference for researchers who wish to perform cost-effective speed congenic screening between a pair of strains from the panel. The assay can be performed in any standard molecular biology labs. Upon publication, this information will be deposited to MGI database that will be available at <http://www.informatics.jax.org>.

ACKNOWLEDGEMENT

The authors acknowledge the review and comments received from Lluís Montoliu. This work was partially supported by an Institutional Development Award (IDeA) to CBG (PI: Shelley Smith) from the National Institute of General Medical Sciences of the National Institutes of Health under grant number P20GM103471 and Center for Humanized Mice Development award from ORIP/DPCPSI/NIH/1R24OD018546-01 (MPI: Gorantla, S/Poluektova, LY). RMQ, DWH and CBG acknowledge the Nebraska Research Initiative and UNMC Vice Chancellor for Research Office for supporting the mouse genome engineering core facility.

References

- Armstrong, N.J., Brodnicki, T.C., Speed, T.P., 2006. Mind the gap: analysis of marker-assisted breeding strategies for inbred mouse strains. *Mamm. Genome Off. J. Int. Mamm. Genome Soc.* 17, 273–287. doi:10.1007/s00335-005-0123-y
- Auerbach, A.B., Norinsky, R., Ho, W., Losos, K., Guo, Q., Chatterjee, S., Joyner, A.L., 2003. Strain-dependent differences in the efficiency of transgenic mouse production. *Transgenic Res.* 12, 59–69.
- Beier, D.R., 2010. New genetic resources for mammalian developmental biologists. *F1000 Biol. Rep.* 2, 72. doi:10.3410/B2-72
- Cao, Y., Liu, X., Deng, N., Jiao, Y., Ma, Y., Hasty, K.A., Stuart, J.M., Gu, W., 2013. Congenic mice provide evidence for a genetic locus that modulates spontaneous arthritis caused by deficiency of IL-1RA. *PLoS One* 8, e68158. doi:10.1371/journal.pone.0068158
- Cendán, C.M., Pujalte, J.M., Portillo-Salido, E., Montoliu, L., Baeyens, J.M., 2005. Formalin-induced pain is reduced in sigma(1) receptor knockout mice. *Eur. J. Pharmacol.* 511, 73–74. doi:10.1016/j.ejphar.2005.01.036
- Collins, S.C., Wallis, R.H., Wallace, K., Bihoreau, M.T., Gauguier, D., 2003. Marker-assisted congenic screening (MACS): a database tool for the efficient production and characterization of congenic lines. *Mamm. Genome Off. J. Int. Mamm. Genome Soc.* 14, 350–356. doi:10.1007/s00335-002-3058-6
- Davie, S.A., Maglione, J.E., Manner, C.K., Young, D., Cardiff, R.D., MacLeod, C.L., Ellies, L.G., 2007. Effects of FVB/NJ and C57Bl/6J strain backgrounds on mammary tumor phenotype in inducible nitric oxide synthase deficient mice. *Transgenic Res.* 16, 193–201. doi:10.1007/s11248-006-9056-9
- Doetschman, T., 2009. Influence of genetic background on genetically engineered mouse phenotypes. *Methods Mol. Biol. Clifton NJ* 530, 423–433. doi:10.1007/978-1-59745-471-1_23
- Estill, S.J., Garcia, J.A., 2000. A marker assisted selection protocol (MASP) to generate C57BL/6J or 129Sv/SvEvTac speed congenic or consomic strains. *Genes. N. Y. N* 2000 28, 164–166.
- Farber, C.R., Corva, P.M., Medrano, J.F., 2006. Genome-wide isolation of growth and obesity QTL using mouse speed congenic strains. *BMC Genomics* 7, 102. doi:10.1186/1471-2164-7-102
- Fleming, R.E., Holden, C.C., Tomatsu, S., Waheed, A., Brunt, E.M., Britton, R.S., Bacon, B.R., Roopenian, D.C., Sly, W.S., 2001. Mouse strain differences determine severity of iron accumulation in Hfe knockout model of hereditary hemochromatosis. *Proc. Natl. Acad. Sci. U. S. A.* 98, 2707–2711. doi:10.1073/pnas.051630898
- Gale, G.D., Yazdi, R.D., Khan, A.H., Lusi, A.J., Davis, R.C., Smith, D.J., 2009. A genome-wide panel of congenic mice reveals widespread epistasis of behavior quantitative trait loci. *Mol. Psychiatry* 14, 631–645. doi:10.1038/mp.2008.4
- Glaser, S., Anastassiadis, K., Stewart, A.F., 2005. Current issues in mouse genome engineering. *Nat. Genet.* 37, 1187–1193. doi:10.1038/ng1668
- Gorham, J.D., Ranson, M.S., Smith, J.C., Gorham, B.J., Muirhead, K.-A., 2012. 1 + 1 = 3: Development and Validation of a SNP-Based Algorithm to Identify Genetic Contributions from Three Distinct Inbred Mouse Strains. *J. Biomol. Tech. JBT* jbt.12–2304–004. doi:10.7171/jbt.12-2304-004

Goto, K., Ebukuro, M., Itoh, T., 2005. Microsatellite-directed selection of breeders for the next backcross generation by using a minimal number of loci. *Comp. Med.* 55, 34–36.

Heiman-Patterson, T.D., Sher, R.B., Blankenhorn, E.A., Alexander, G., Deitch, J.S., Kunst, C.B., Maragakis, N., Cox, G., 2011. Effect of genetic background on phenotype variability in transgenic mouse models of amyotrophic lateral sclerosis: a window of opportunity in the search for genetic modifiers. *Amyotroph. Lateral Scler. Off. Publ. World Fed. Neurol. Res. Group Mot. Neuron Dis.* 12, 79–86. doi:10.3109/17482968.2010.550626

Iakoubova, O.A., Olsson, C.L., Dains, K.M., Choi, J., Kalcheva, I., Bentley, L.G., Cunanan, M., Hillman, D., Louie, J., Machrus, M., West, D.B., 2000. Microsatellite marker panels for use in high-throughput genotyping of mouse crosses. *Physiol. Genomics* 3, 145–148.

Inoue, M., Sakuraba, Y., Motegi, H., Kubota, N., Toki, H., Matsui, J., Toyoda, Y., Miwa, I., Terauchi, Y., Kadowaki, T., Shigeyama, Y., Kasuga, M., Adachi, T., Fujimoto, N., Matsumoto, R., Tsuchihashi, K., Kagami, T., Inoue, A., Kaneda, H., Ishijima, J., Masuya, H., Suzuki, T., Wakana, S., Gondo, Y., Minowa, O., Shiroishi, T., Noda, T., 2004. A series of maturity onset diabetes of the young, type 2 (MODY2) mouse models generated by a large-scale ENU mutagenesis program. *Hum. Mol. Genet.* 13, 1147–1157. doi:10.1093/hmg/ddh133

Johnson, K.R., Zheng, Q.Y., Noben-Trauth, K., 2006. Strain background effects and genetic modifiers of hearing in mice. *Brain Res.* 1091, 79–88. doi:10.1016/j.brainres.2006.02.021

Lamacchia, C., Palmer, G., Gabay, C., 2007. Discrimination of C57BL/6J Rj and 129S2/SvPasCrl inbred mouse strains by use of simple sequence length polymorphisms. *J. Am. Assoc. Lab. Anim. Sci. JAALAS* 46, 21–24.

MacLennan, M.B., Anderson, B.M., Ma, D.W.L., 2011. Differential Mammary Gland Development in FVB and C57Bl/6 Mice: Implications for Breast Cancer Research. *Nutrients* 3, 929–936. doi:10.3390/nu3110929

Markel, P., Shu, P., Ebeling, C., Carlson, G.A., Nagle, D.L., Smutko, J.S., Moore, K.J., 1997. Theoretical and empirical issues for marker-assisted breeding of congenic mouse strains. *Nat. Genet.* 17, 280–284. doi:10.1038/ng1197-280

Masuya, H., Shimizu, K., Sezutsu, H., Sakuraba, Y., Nagano, J., Shimizu, A., Fujimoto, N., Kawai, A., Miura, I., Kaneda, H., Kobayashi, K., Ishijima, J., Maeda, T., Gondo, Y., Noda, T., Wakana, S., Shiroishi, T., 2005. Enamelin (Enam) is essential for amelogenesis: ENU-induced mouse mutants as models for different clinical subtypes of human amelogenesis imperfecta (AI). *Hum. Mol. Genet.* 14, 575–583. doi:10.1093/hmg/ddi054

Montagutelli, X., 2000. Effect of the genetic background on the phenotype of mouse mutations. *J. Am. Soc. Nephrol. JASN* 11 Suppl 16, S101–105.

Ogonuki, N., Inoue, K., Hirose, M., Miura, I., Mochida, K., Sato, T., Mise, N., Mekada, K., Yoshiki, A., Abe, K., Kurihara, H., Wakana, S., Ogura, A., 2009. A high-speed congenic strategy using first-wave male germ cells. *PloS One* 4, e4943. doi:10.1371/journal.pone.0004943

Petkov, P.M., 2004. An Efficient SNP System for Mouse Genome Scanning and Elucidating Strain Relationships. *Genome Res.* 14, 1806–1811. doi:10.1101/gr.2825804

Petkov, P.M., Cassell, M.A., Sargent, E.E., Donnelly, C.J., Robinson, P., Crew, V., Asquith, S., Haar, R.V., Wiles, M.V., 2004. Development of a SNP genotyping panel for genetic monitoring of the laboratory mouse. *Genomics* 83, 902–911. doi:10.1016/j.ygeno.2003.11.007

Puccini, J., Dorstyn, L., Kumar, S., 2013. Genetic background and tumour susceptibility in mouse models. *Cell Death Differ.* 20, 964–964. doi:10.1038/cdd.2013.35

Sakai, T., Miura, I., Yamada-Ishibashi, S., Wakita, Y., Kohara, Y., Yamazaki, Y., Inoue, T., Kominami, R., Moriwaki, K., Shiroishi, T., Yonekawa, H., Kikkawa, Y., 2004. Update of mouse microsatellite database of Japan (MMDBJ). *Exp. Anim. Jpn. Assoc. Lab. Anim. Sci.* 53, 151–154.

Siuciak, J.A., McCarthy, S.A., Chapin, D.S., Martin, A.N., Harms, J.F., Schmidt, C.J., 2008. Behavioral characterization of mice deficient in the phosphodiesterase-10A (PDE10A) enzyme on a C57/Bl6N congenic background. *Neuropharmacology* 54, 417–427. doi:10.1016/j.neuropharm.2007.10.009

Suemizu, H., Yagihashi, C., Mizushima, T., Ogura, T., Etoh, T., Kawai, K., Ito, M., 2008. Establishing EGFP congenic mice in a NOD/Shi-scid IL2Rg(null) (NOG) genetic background using a marker-assisted selection protocol (MASP). *Exp. Anim. Jpn. Assoc. Lab. Anim. Sci.* 57, 471–477.

Teppner, I., Aigner, B., Schreiner, E., Müller, M., Windisch, M., 2004. Polymorphic microsatellite markers in the outbred CFW and ICR stocks for the generation of speed congenic mice on C57BL/6 background. *Lab. Anim.* 38, 406–412.

Thifault, S., Lalonde, R., Sanon, N., Hamet, P., 2002. Comparisons between C57BL/6J and A/J mice in motor activity and coordination, hole-poking, and spatial learning. *Brain Res. Bull.* 58, 213–218.

Tsang, S., Sun, Z., Luke, B., Stewart, C., Lum, N., Gregory, M., Wu, X., Subleski, M., Jenkins, N.A., Copeland, N.G., Munroe, D.J., 2005. A comprehensive SNP-based genetic analysis of inbred mouse strains. *Mamm. Genome Off. J. Int. Mamm. Genome Soc.* 16, 476–480. doi:10.1007/s00335-005-0001-7

Wakeland, E., Morel, L., Achey, K., Yui, M., Longmate, J., 1997. Speed congenics: a classic technique in the fast lane (relatively speaking). *Immunol. Today* 18, 472–477.

Witmer, P.D., Doheny, K.F., Adams, M.K., Boehm, C.D., Dizon, J.S., Goldstein, J.L., Templeton, T.M., Wheaton, A.M., Dong, P.N., Pugh, E.W., Nussbaum, R.L., Hunter, K., Kelmenson, J.A., Rowe, L.B., Brownstein, M.J., 2003. The development of a highly informative mouse Simple Sequence Length Polymorphism (SSLP) marker set and construction of a mouse family tree using parsimony analysis. *Genome Res.* 13, 485–491. doi:10.1101/gr.717903

Table 1. Primer sequences used for PCR amplification of SSLP markers

Marker	cM	Primer F	Primer R	Product	Lab
				Size	Code
D1mit231	6.5	ACCCACAATTGCCTGTGG	GTCTTTGCAAGCCACCAAAT	267	1-1
D1Mit211	10.59	GTTATTCATCAAAATACAGATGGCC	TCTGCTGCTAAGTAGAATGAATGC	135	1-2B
D1mit213	22.88	TTCTTAGAAGTGATAAAAGTTTCAGCA	AAAATTCCAGAATTCTCACTACGG	108	1-2
D1mit303	31.79	GGTTTCTATTTCCGTTCTCGG	TCTGTGCTGCAAAACAGAGG	128	1-3
D1mit46	39.16	AGTCAGTCAGGGCTACATGATG	CACGGGTGCTCTATTTGGAA	253	1-4
D1Mit440	44.98	TCCACACAAGGTGTCCTCTG	GCTCAGGTGACCTCCAAAAC	114	1-5B
D1Mit218	55.76	TGCAAAATGTACTTTAGTCTCTAGTGG	AGTTTTGGTGAGTGATCTTATCCC	145	1-218
D1Mit369	63.51	ACTTGTGTTGTTGCTGAGGTTCA	GCTTATGAACCCACCCTCAA	141	1-5Q
D1Mit200	63.84	GCCATGTTTATGTACATAGGTAGG	ATGGATGGATGGTTTTCTG	199	1-200
D1mit507	72.38	GTTGGAAAGACTGTTTACAATTCTG	GTTCCAGCCTTTGCCCTTAC	115	1-6
D1Mit150 F	81.08	CTGGTGTCACACACAGGTC	TGATTGCAATATACCAGGTTTCC	138	1-150
D1Mit407 F	89.89	GAGAACAACCGCCACCAAT	ATATTTGCTTTGAAGTTACTTTGTGTG	119	1-407
D1Mit155	98.2	ATGCATGCATGCACACGT	ACCGTGAAATGTTCACCCAT	252	1-7
D2mit1	2.23	CTTTTTCGTATGTGGTGGGG	AACATTGGGCCTCTATGCAC	123	2-1
D2Mit149	10.04	ATATCATATAGTAGAGAAAGCGTGCTG	TCATTAGACTTGGAAAAAGTTTGC	199	2-1P
D2mit152	21.81	CACAGATCTTGTAAAGACCACGTG	TGCCATGAGTGTGGGACTAA	109	2-2
D2Mit367	22.5	GCCTGTGCTAAAAAGAGGTG	GCCCTGAGAACTACCCTCCT	149	2-367
D2Mit323	31.42	AGAATCCTAAGTGGTGGTTAGAGG	ACCCAAAGTTGCTTTAAGTACACA	125	2-2P
D2mit328	42.89	CTTTCAATGTTCCGGCATG	AAGACTTGCTTTCATTAGACCACA	235	2-3
D2Mit42	54.85	ATTACTGGGCAGGAACATTTG	GCCAACTTCCAGACTCCTC	132	2-3B
D2mit395	59.97	AGGTGAGCCTGGACTATATGG	AGCATCCATGGGATAATGGT	125	2-4
D2Mit106	64.78	GAGGGTTGCCAAGAGACTG	CACCTCAGGGGAACATTGTG	150	2-4P
D2Mit107	65.13	GGGAGTGAAGCCAGCATAAG	AACTGACTGAGTTTCAAAGTGCC	119	2-107
D2Mit309	75.03	ACAAATGCCACTCTCACATCC	TATTTCTCAGAGTCACTAGGAGTGATG	119	2-309
D2mit285	75.41	TCAATCCCTGTCTGTGGTAGG	TATGACACTTACAAGTTTTTGGTG	141	2-5
D2Mit48	77.36	GCTCTGCAGAAGATGCTGC	GCTGAGACGCAGAGTCGC	130	2-5B
D2Mit311	86.02	ACAGGCAGCCTTCCCTTC	TCTGTCCCCTCTGTTTCT	126	2-5D
D2Mit265	97.97	AATAATAATCAAGTTGTCTTGAACC	TAGTCAAAATCTTTTGTGTGTTGC	105	2-6P
D2mit266	103.83	GGATCTATGCTCATTTTAATTGC	TCATCTTCTGGTTTCAACATGG	127	2-6
D2Mit148	100.49	GTTCTCTGATCTACGGGCATG	TTCACTTCTACAAGTTCTACAAGTTCC	115	2-148
D3mit60	1.96	GACATCCTGGGCAACATTG	GGTGTTGTTTGCTGTTGCTG	170	3-1
D3mit164	2.01	GCTCCTGGGAAAGGAAGAAT	GATACTTGGGGTTGTGCATACA	135	3-164
D3mit203	10.82	CTGAATCCTTATGTCCACTGAGG	GGGCACCTGCATTCATGT	150	3-2
D3mit51	35.2	GGCACTGATAGCAGGCCTAG	TCTCTTCTGGTATTTCTCTCCG	248	3-3
D3mit 49	39.02	CTTTTCTGCCCCACTTTC	TCCTTTTAGTTTTGATCCTCTGG	132	3-49
D3Mit10 F	43.03	CTGGCTTGGTGAAGTCTT	CCTAAGCCAGCTACCACCAC	142	3 10
D3Mit189	43.89	GTTACCACCAGAGAAAAGGC	TACTCTCGCTGCTTCCCTA	133	3-189
D3Mit318	54.58	CTCATTCCTTCTGAGCAATGG	TATGGGATATGCTTTTCATAAAAGG	148	3-318
D3mit256	63.04	TACATTGCTTTTTGCTTTGAGTG	GTCGAATGTTATCAGAATTTGCA	125	3-4
D3mit116	79.23	TCACTGCCCATCTTTGTAACC	CCCAGAGACCCGGAATAGAA	259	3-116
D3mit19	87.6	CAGCCAGAGAGGAGCTGTCT	GAACATTGGGGTGTTTGCTT	159	3-5

Table 1. Continued: Primer sequences used for PCR amplification of SSLP markers

Marker	cM	Primer F	Primer R	Product Size	Lab Code
D4mit227	4.43	CTCAGACATGATTTTTTCCAAGG	GCAGTTAAACTGTACTTTCTGTAAACA	181	4-1
D4Mit193	13.99	TATTTTAATTTTAGCCCATCAGGG	AAAGACATACAATTGATCCACAGG	136	4-1B
D4mit17	33.96	TGGCCAACCTCTGTGCTTCC	ACAGTTGTCCTCTGACATCC	147	4-2
D4Mit27	42.13	GCACGGTAGTTTTTCCAGGA	TGGTGGGCAGGCAATAGT	150	4-2B
D4Mit31	50.04	ACGAGTTGTCTCTGATCAACA	AGCCAGAGCAAACACCAACT	121	4-2C
D4mit308	57.66	TATGGATCCACTCTCCAGAAA	CAAAGTCTCTCCAAGGCTG	88	4-3
D4Mit203	63.26	GAATTCTTCCTGGGCTTTC	CAAGAGCCAGGTGTGGTAT	105	4-3B
D4mit251	69.05	AAAAATCGTCTTTGACTTCTACATG	TTTAAAGGGTTCTTTATCCTGTG	114	4-4
D4Mit354	76.12	TTGATCTGTCTGGATTCCA	AGACAGACACATAGACAGACATAGG	111	4-4B
D4mit42	82.64	CATGTTTGCCACCCTGAAAC	CCTCACTTAGGCAGGTGACTC	100	4-5
D5mit145	3.37	TATCAGCAATACAGACTCAGTAGGC	TGCCCCCTTAAATTCATGGTC	150	5-1
D5Mit73	12.88	GTTTGAGAGGTCTGAAAGCA	TTTCCATTTACATACATTTGTGCA	113	5-1B
D5mit352	18.4	CCCAGAGCCCACATCAAG	TAGGTGGGTGTGTCTCTCCC	110	5-2
D5Mit233	28.55	TCCCCTCTGATCTCCTCAGA	CCTCTAGAATACAATTCAATGTGG	147	5-3
D5Mit307	40.31	AGTGGCATAAATTCATTCAATAATG	GGAATCAAGTTGTTTTTAAATTTACC	120	5-3P
D5Mit277	51.7	GTGTGTTTGTGCATGGGTATG	ACCATCGGGAAAAAATGTAGC	123	5-4B
D5Mit161	65.34	CACACGCATAGTCTTGTTGGG	GCATGTTCAACTGTGCTTTCA	119	5-5B
D5Mit99	79.51	CAGAAAAGAGAAAACGGAGGG	TTCTGCTGCCTGAAGTTT	200	5-5C
D6mit138	1.81	GCTCTTATTAATGAAGAAGAAGGAGG	CAAAGAAAGCATTTCAGACTGC	135	6-1
D6Mit296	2.25	TCGGGCATCTTTATTTTTGC	TAGTGACGACACCCCCCT	100	6-296
D6Mit116	11.5	ACATTTCTTTGTGAGGTTCCCTTG	CAGGTTTTTTGAAAGACACTCTTG	122	6-1B
D6Mit274	23.7	GCAATGCCAAAATGTTCAAA	TCCTTCTCCATTTACACTTACAACA	113	6-1D
D6Mit8	35.94	TGCACAGCAGCTCATTCTCT	GGAAGGAAGGAGTGGGGTAG	163	6-2B
D6mit100	41.03	CTTGAGTAGGTCTCAGTGCGG	CACATGCACACACAGAAGCA	84	6-3
D6mit36	48.93	ACCATCTGCATGGACTCACA	GTTGAAGAGGACGACCAAGTG	196	6-4
D6Mit10	52.75	TCAGAGGAACAAAGCAGCAT	CCTGTGGCTAACAGGTAAAA	200	6 10
D6Mit333	60.55	TCCTCACTACAATTCATCTATTACTGC	TGCTTCTGGTATAGGCAGTTAGG	122	6-333
D6Mit302	67.66	AATGACCCTGGTTAGTGTGAGG	GAATTCCATTGAGGGGC	113	6-4C
D6mit14	77.64	ATGCAGAAACATGAGTGGGG	CACAAGGCCTGATGACCTCT	157	6-5
D7Mit21	1.91	GGGTGAACCTTACAGGGGT	ATCAAACCAGCCCAAGTGAC	192	7-0A
D7Mit152	2.71	GCCTAGCACACGCCAAAG	CCTTGTGCATGTTGCTATG	129	7-152
D7mit57	9.94	TTCCCTCTAGAACTCTGACCTCC	AGTTCAGAGCCGAGACTAGGC	148	7-1
D7mit267	17.09	CTCTTTCTGTTACATGGTTAGATTCC	AAAGACAGTTGAAGTTGACTTCTGG	192	7-2
D7Mit158 F	29.82	CTTCATCTGAGCCTGGGAAG	ACTGTAGACCCATGTTCTGATTAGG	149	7-158
D7Mit317	41.5	ATGTCTCCTTGACATTGGGC	TCTTGTAATCTCACATCTAAGTGTGTG	102	7-2D
D7Mit194	44.83	GTGCAACACACAGAAAAGTTCC	AAAGGCTCACAACGGACTGT	146	7-194
D7Mit220	55.69	AAGCATGCAAGCACACTCAC	ATGCACACAGGCAGTCACTC	135	7-3A
D7mit101	69.01	TACAGTGTGAACATGTAGGGGTG	TCCCAACATGGATGTGCTAA	111	7-4
D7mit223	88.85	ATGCACATGAGTGTGTGTATGC	TCCTGTGTCTGACGCTCATC	106	7-5

Table 1. Continued: Primer sequences used for PCR amplification of SSLP markers

Marker	cM	Primer F	Primer R	Size	Code
D8Mit143	11.59	AGCCTGAGGTTATGTTTTTGC	GGCCCTCAGGTTTCTCTCT	279	8-143
D8mit178	34.43	AAAATCAACTGTTTACATTTGAGCC	AGAGCACGCAGTGTGTATGC	148	8-2
D8mit45	42.16	GAACAGGACCAATAAAATGAAAGC	CTACCTTACCAAACCTCCCGG	121	8-45
D8Mit80	43.06	TGCATTTGTCAGGGCTCTC	ATGACACATGAGCCTCCACA	107	8-080
D8mit211	52	CAGAACTGTCCTGAAAAGTCC	TACCCACAAACCTGTATTTAAATTAA	149	8-4
D8mit215	62.63	AATACACAAGTTGGCCTCA	ATGTGTGGATATTCATGTGCTC	178	8-5
D8Mit56 F	76.14	ACACTCAGAGACCATGAGTACACC	GAGTTCACTACCCACAAGTCTCC	162	8-56
D9mit250	2.46	CCCCAAAACCTATTTGCACTG	GTGACATGATTCTTCAGTCTTACC	123	9-1
D9Mit249	2.46	AAGCCCTCTTAGAAGTAGTGTATG	AGCCATGAACTAATTACATGTATCA	125	9-249
D9Mit89 F	14.79	CACATACAAGGATATACATACAGGC	TCACAGGAGGTGGCAGAAAT	145	9-89
D9mit285	21.42	CAATACATTGCTGATTATATCAGAGA	GGACTCTAGATCTCATCAGGGA	125	9-2
D9mit336	35.39	AAGTGGTTCACAGAAATGTATACAGG	TTTTCTTTCTGTGGTAAAGGGG	122	9-3
D9Mit8	42.65	GATGAAGACAATAAAGAACCTTAAA	AAGAGCTAACCCATTGCTGC	183	9-3B
D9mit347	55.11	CCTCCACATGTGCACTGCT	CTGTCCATCTATCATCTATCTGTCTG	122	9-4
D9mit18	71.49	TCACTGTAGCCAGAGCAGT	CCTGTTGTCAACACCTGATG	180	9-5
D10Mit28	3.04	CCTCCTGTATGTGATTTAAAGCA	CTGCCATCTGACCCTGATA	147	10-1A
D10mit213	9.75	CTCCTCTACTGATTGTCCCC	GGGACAACTTTTAAAAATTGCA	150	10-2
D10Mit108	22.89	TGCCTGTAACCTGCATACCC	GTTTAACACCCAGGACTATACATGG	142	10-108
D10mit20	34.83	CACCCTCACACAGATATGCG	GCATTGGGAAGTCCATGAGT	234	10-3
D10mit230	45.28	AGATAGCCTAGGGGGTGCAT	ATCAGTTTCCAATCGTGCT	115	10-4
D10Mit96	51.19	CTTCTTTGAAGTTAGATGCAGCC	TACGGAGAAGGGAACACCTG	150	10-096
D10Mit178	51.42	ATTGTCAAATATCTTCTCAGTTGC	TTATTCCTAGGCAGTCTGTCTGG	133	10-4B
D10mit233	61.58	GTGCTTTATATTGGAGATCATACA	GTCCCGAATTTACATACATAGC	130	10-5
D10mit271	72.31	ACAACCAAAGTCTTTGTAGAAGA	AATATATAGGCACACCTTAATAGCCA	117	10-6
D11Mit1	2.51	GGGTCTCTGAAGGCTTTGTG	TGAATACAGAAGCCACGGTG	153	11-1P
D11mit226	5.64	AGGTGAACTCTTTGAAGTTTGTG	AAAGGAGTGACTGAGAAAGACACC	139	11-226
D11Mit151	15.29	TGGGAATTCTGGGAGTTCTG	GTTGGTCTGTTATGAAGACCAGG	140	11-1Q
D11mit217	23.1	ACTGGAAAATATGTTTAAACCTCTG	AAATGGGATTCTGCAAAAACC	135	11-217
D11Mit349	33.29	AGTATCAGAAGATCCAGTTGGAGG	GTAGAAAAAGATACCCAGTGTGAGC	118	11-349
D11Mit298	42.76	AAACAAACAAAAATGCACCTCA	GTACCACCATGCCTAGCCTC	199	11-298
D11mit39	52.22	TTTCATGACCCCTAATTTCCC	GTGGGTGTGCCTGTCAATC	155	11-39
D11mit70	58.9	GGAAGTAGCTATGGAGGTGGC	TCTGACCCAGAGCTCAAATACA	140	11-70
D11mit54	59.82	AGGCTGGTGGCTAGTGTCC	AAGTCTTGGCTGCATCTTT	144	11-4A
D11Mit132	62.92	GGTCAGAGGACAATCTTACATGC	GTTCCAAGACAATGAGAGACCC	117	11-4B
D11mit333	71.83	CATGTGGTTATTTTCTAGCCCC	AGGCATCAATAACTATTTTTCAGTG	125	11-5
D11mit303	82.9	TCAATCTCTCAAGTTTTTCCAGC	GACAACGTTGACCTCCACG	106	11-6
D12mit12	8.49	TTCAATGCCTTCTGGCTTCT	GATTACCGGGTGTGTGACCT	145	12-1
D12mit2	18.94	ACACAGGCTAAACATGGGC	GCATCTGTATTCCACAGGCA	134	12-2
D12mit114	28.94	TTGACCTTGAACCTGTGACCC	GTTTTCTCCAAATCACTGTCACC	144	12-3
D12Mit149	36.85	CATGGCACACATACATACATGTG	AACATAGCAATGGTATATAGGTATGGG	132	12-3B
D12Mit117	43.32	AATTGAGGAACTTAGAAGAAAAGCC	CCTCTGGCCACCATAATG	127	12-117
D12Mit30	50.43	TATGTGACTGCAATCCCAGC	ATGAACACATCATGCCAGA	107	12-30
D12Mit99	52.9	CTTACAGAAAATGAAAACCAAAACA	CCTCTGCTTTAGAGGCAAACG	151	12-099
D12mit263	62.11	TCAGATCTCAGCAGATAAATACTTGG	TCCCTGGAGCATATTTGAC	113	12-6
D12nDs2	62.22	ACATGGTAATTTATGGGCAA	CTGGATACCTGCAATAGTAGA	195	12-5

Table 1. Continued: Primer sequences used for PCR amplification of SSLP markers

Marker	cm	Primer F	Primer R	Size	Code
D13mit16	7.26	CCAGCTGAAGGCTTACTCGT	AAAGTTAGAATCAGCCATTCAAGG	207	13-1
D13Mit137	18.87	GAATCAGAGAACTGGCTGTG	TCTAAAAGAGAGAAAATTGGGGC	131	13-137
D13Mit63	21	GAGATGGAAGGAAGAGATGGC	CAAATGCATGTATCCGTATGTG	139	13-63
D13mit186	31.87	GAAAGCCCTAGGGGAAGATG	TGCAGTTTCTAAGGTTAAACTAAAGC	149	13-3
D13Mit282	32.59	TCGCACTTCTATACAGTTATAAGAG	GGACAGAAAGCATGCAGAGG	125	13-282
D13Mit191	45.05	GCAAATTAGAGAATCATGCATCC	TGGAGAATGCTAAAAGCATGG	119	13-3B
D13Mit51	56.45	TCCTGCAAAAGTGGAGCC	TGGAAACAAGCTCTTGGAGG	146	13-3P
D13mit213	59.69	GCCTGAAACTCTACATAAAATACATCC	AGTTTCATTGCTTTAGTTACATTTTCA	149	13-5
D13mit78	67.21	ACAGCACGGGTTTATCATCC	TATGCCTGCCAGGCTTCTAT	229	13-6
D14Mit179	4.92	CCACTTGCAGCATTGACAAT	CAAACATCTGTGACAATAAAATTTCA	144	14-1P
D14mit126	11.94	CCTGTCCACAACACCTTTT	TATACATATGGGTAGCACTGAGTGG	140	14-2
D14mit60	24.6	AGGCTGCCCATAAAAGGG	GTTTGTGCTAATGTTCTCATCTGG	132	14-3
D14Mit259	27.65	TGGTGTCTCCTTCGGAATTT	TAAATGTAAGGTAAAGGCAATGG	125	14-259
D14Mit39	35.69	AAAGAGCAACCCCAATTCT	ACTTTTACCTGGTCTCCAAAAGC	246	14-39
D14mit68	37.61	GTGGCATGCACAACCGTATA	CCCTTTTGAGGTGCTTGT	153	14-4
D14Mit194	45.96	AATATTCTAAATGAAATCCAATGTGTG	TTAATTGCAAGTAACACAATGAGTAGG	92	14-4B
D14mit95	57.2	TATTTTTAAGTCAGTATACACATGCGC	TTATCCAAGTGATTTAAAGAAGAGGC	124	14-5
D14Mit97	62.2	TCAGTCCAACTCTGTTAATCTTCC	CAGCTCCACATTTTGTCTCA	136	14-6A
D15Mit13	1.84	GGAGACAAAATGAACCTCTGG	TTGTAAGACAAGCATAGCTCAACA	136	15-13
D15Mit265	6.08	ACATTAGTCAACTATGCTGGTACTCTG	TTCTCTCTAATGTCAATTGTTTCA	198	15-265
D15Mit138	15.68	TTCAATTCCTTTTGTCAAATG	CAAGACCTTAGATTCACTCTACCC	147	15-138
D15mit154	16.82	AGCACTGGGTACACAACTGG	ATGAAAGCATGTGTAGTCTTTCTCA	150	15-2
D15mit128	25.58	CAAGTTCTGCAAGAATTATTATGC	CCACTTTGAAGTTTTCTTTCTTAGC	134	15-3
D15Mit92	32.19	AGTCTCTCTCCCCCTTCTCTC	TGCCACAAGCACAAATAGTATCC	147	15-3B
D15mit80	38.02	CATTGAGGGTTGTAGGTTGG	ACCCCTGCAAGTTGTCTTTG	149	15-4
D15Mit34	45.31	TGGACAACCATTTTGGACAA	CTTTCTGTGAGGCATCACCA	147	15-34
D15Mit244	48.65	TCTACCCTCTGTGGAACATCG	CTTTGTGTCCATACACTAATATCAAGG	116	15-4P
D15mit16	58.05	AGACTCAGAGGGCAAAATAAAGC	TCGGCTTTTGTCTGTCTGTC	119	15-6
D16mit131	3.41	TGGTGGTGGTGTGATGGTA	AAGACCATTCTAATAAACAACACCC	140	16-1
D16mit136	27.82	AGATAATTCCCGTGAGAATAAAACC	TTGAGAAGTTTGCCTATAATGG	129	16-2
D16Mit30	33.01	GTGCACATACATACCACAGCG	TCACTGCAGGGAGGTTTCTAG	152	16-30
D16mit64	34.22	TACCATGATCAGTCCAAAGGC	ACTTAAGGTTGTCTGTGGGG	220	16-3
D16Mit140	40.3	ATAGTTGAAAACTTGAACATGCG	GAAAAGGTTAATGCTGGTCACC	145	16-3B
D16Mit19	45.36	CAGGCATGTGAACAAAGTGG	GTGACTGATGAATGCCTGACA	121	16-3C
D16mit70	48.81	GGATCTATATGCTATAGAACCATTCA	GTCATCAATTCATTTCTAATATAGA	187	16-4
D16Mit71	57.06	TAGAAAATCTTCAAATAGGATCTGTTT	GAGCATTTCCCTTTTACCTGG	154	16-6A
D17mit164	2.11	AGGCCCTAACATGTAGCAGG	TATTATTGAGACTGTGGTTGTTGTTG	133	17-1
D17Mit133	12.53	TCTGCTGTGTTACAGGTGA	GCCCCTGCTAGATCTGACAG	188	17-1C
D17mit51	19.74	TCTGCCCTGTAACAGGAGCT	CTTCTGGAATCAGAGGATCCC	154	17-2
D17Mit68	23.55	GTCCTGACATCATGCTTTGTG	CTACCGTTTGGAAGGCTGAG	130	17-68
D17mit20	29.73	AGAACAGGACACCGGACATC	TCATAAGTAGGCACCAATGC	165	17-3
D17mit205	39.3	TGTGCATGTATATGTGTGTAATG	GCTAAGTCAGAAGAGTTCTGTAATGG	223	17-4
D17Mit1002	50.97	TCTGAATGCTGACTTCCATCC	ACACATATGTGTAGTGTATGAATGTGC	138	17-4P
D17mit123	60.67	CACAAGGAGGGAGCCTGTAG	CACCGTAAGAGTCTAATAATAAGGGG	133	17-6

Table 1. Continued: Primer sequences used for PCR amplification of SSLP markers

Marker	cM	Primer F	Primer R	Product Size	Lab Code
D18mit222	8.08	AATCCAAGATTGACATGTGGC	CTTAGATGCCCTGTCTTAAAAAA	113	18-1
D18Mit226	18.18	CAGGCAGGGTGCATATATTATAA	TATCTGTTTATGTGTGTACATTGTGTG	125	18-1C
D18mit177	21.39	CTGTAGTTTATCAGTTCACCCTGTG	TGTGCTGTAAACAAATATCTCTGG	171	18-2
D18mit91	29.67	TCCACAAATGTTGGCAAAGA	TTTCTGGCCATATTGGAAGC	140	18-3
D18mit124	32.15	CCCAATGGGGTGTCTTTTA	CTGCCACACATTGTGTGTATG	150	18-4
D18Mit184	39.7	CACACATGTGTAGGTAGGTAGGTAGG	CGCACAAAGGACTACTGAAACA	172	18-4C
D18mit186	45.63	AAGTGTGGGCAAAGGCTAA	CTTTAGTATAGTGTGCATGAGTGTGA	125	18-5
D18Mit7	51.92	ACAGGAGAACGGGAAGCTAG	GCCAGAGTGGACCAAGATGA	95	18-5B
D18Mit129	53.28	CCAGCACAGAGGCAGTCAT	TGATTCTTGGGTCTCTGAATACA	138	18-5C
D19mit68	3.38	CCAATACAAATCAGACTCAATAGTCG	AGGGTCTCCCATCTTCTTA	132	19-1
D19Mit60	13.9	CAACACCTCACTGTTTAGGAACC	GCTGAGGTCAATATTTAGCATGC	136	19-60
D19Mit45	16.14	CCATTCATAAAATGGGCTTAGG	ACCATGAATGTGTTTTGAGGTG	145	19-2A
D19Mit30	21.34	GGTGGCTTAGAAATAGTATCGAAA	CCAGCTCTAGGCAGGCATAT	150	19-2B
D19mit88	32.23	AACAGTGCAACTTTGGAGGC	TCATTGGAAGTGTCTTAACAGTGC	148	19-3
D19Mit19	34.08	CCTGTGTCCATACAGGCTCA	ACCATATCAGGAAGCACCATG	138	19-19
D19mit11	36.26	TCAAAGTCAAGGTGGGCAG	ACTTTCAGATGTTGGGCAC	146	19-4
D19mit1	50.32	AATCCTTGTTCACTCTATCAAGGC	CATGAAGAGTCCAGTAGAAACCTC	122	19-5
D19mit71	56.28	ATGATTCCCGCAGTTTGTGTC	TCTCAACTGTTATTCCTCAATAGCC	136	19-6
DXMit136	4.23	ACGGAAACACTCTTATGTGCG	ATTTTGATTACAGCATGTCCCC	182	X-136
DXMit48	25.51	ACCTCCCCCTGCATTACTCT	TTCTCCAGAATCCATGCTCC	105	X-48
DXMit62	34.6	GCAATTGTGATGTTGTAGTAAATATGG	ATAACTGAGGTCTGCGGGG	125	X-62
DXMit110	35.53	TGACATGAAGTATGTGTTCTGTC	TAGGCACATGTTACATGGG	135	X-110
DXMit170	45.87	TGCAGGCACTAACAGTGAGG	TAGTTTCACTGTGCCATTGTATACA	115	X-170
DXMit179	53.17	TTTGATAGAGCCATGTTTGGC	CAGGCTAGCCTCAAACCTCTCA	122	73.26
DXMit79	68.46	AGTCTGCCTTCTCTTTCTGTATCC	TGAAACTATTCCAACATTATTCTTGG	138	53.17
DXMit153	73.26	CAATCAAGCAGATGGAAGCA	AAGGACTGCCAAGAGGACAA	143	68.46

Figure 1. Agarose Gel images of PCR bands and their interpretation. Bands were assigned numbers 1, 2, 3 or 4 to indicate their sizes compared to the rest of the bands in that set. Number 1 was assigned to the smallest sized polymorphic band. 2, 3 and 4 represent the successive bigger sized bands.

Band size Interpretation

Marker	cM	Gel Image	Band size Interpretation									
D1Mit231	6.5		1	2	2	2	1	1	1	1	1	2
D1Mit211	10.59		2	1	2	2	2	2	1	3	2	2
D1Mit213	22.88		1	2	1	1	1	1	1	2	1	1
D1Mit303	31.79		1	2	1	1	1	1	1	2	1	1
D1Mit46	39.16		3	1	2	2	2	2	2	1	3	3
D1Mit440	44.98		2	1	2	2	1	1	1	1	1	2
D1Mit218	55.76		2	1	2	3	2	1	1	2	2	1
D1Mit200	63.84		3	1	1	1	2	3	3	1	3	1
D1Mit369	63.51		2	1	1	1	1	1	3	2	2	3
D1Mit507	72.38		1	3	3	3	2	3	4	3	3	4
D1Mit150	81.08		2	3	1	1	1	1	2	1	1	1
D1Mit407	89.89		2	1	2	2	2	2	2	2	2	2
D1Mit155	98.2		1	2	1	2	1	1	1	1	1	1
D2Mit1	2.23		3	1	1	1	1	1	2	1	1	1
D2Mit149	10.04		3	3	1	2	3	2	3	2	2	2
D2Mit152	21.81		2	1	2	2	1	1	1	1	1	1
D2Mit367	22.5		2	1	1	1	2	1	1	1	1	2
D2Mit323	31.42		3	3	1	2	3	2	3	2	2	2
D2Mit328	42.89		3	2	1	1	2	2	2	3	2	3
D2Mit42	54.85		1	3	1	1	4	2	1	4	3	1
D2Mit395	59.97		3	1	3	3	2	1	3	2	1	1
D2Mit106	64.78		1	1	2	2	1	2	2	2	2	2
D2Mit107	65.13		1	1	1	1	1	2	1	1	2	1
D2Mit309	75.03		1	1	3	3	1	1	1	1	1	2
D2Mit285	75.41		3	1	3	3	2	2	2	2	2	1
D2Mit48	77.36		1	2	1	1	2	1	1	1	1	1
D2Mit311	86.02		1	4	3	3	1	2	1	2	1	2
D2Mit265	97.97		3	1	2	2	2	3	3	3	2	1
D2Mit148	100.49		1	2	3	3	2	3	4	3	2	2
D2Mit266	103.83		1	2	1	1	1	1	1	1	1	1

Figure 1. Continued.
Agarose Gel images of PCR bands and their interpretation.

Marker	cM	Gel Image	J	J	J	2	J	J	J	C	e	D
D3Mit60	1.96		1	2	1	1	1	1	1	2	2	1
D3Mit164	2.1		1	2	2	2	1	1	1	1	1	1
D3Mit203	10.82		1	2	1	1	2	1	2	3	1	2
D3Mit 49	21.6		1	2	2	2	3	1	1	3	1	3
D3Mit51	35.2		3	2	2	2	3	1	2	2	1	2
D3Mit10	43.03		3	3	2	2	3	2	3	2	2	1
D3Mit189	43.89		1	1	1	1	2	2	1	2	2	2
D3Mit318	54.58		1	1	2	2	1	3	1	3	3	1
D3Mit256	63.04		1	2	3	3	1	1	1	1	1	2
D3Mit19	87.6		2	1	1	1	3	3	2	3	3	3
D3Mit116	79.23		2	1	1	1	2	1	2	2	2	2
D4Mit227	4.43		1	2	1	1	2	2	3	2	3	3
D4Mit193	13.99		1	2	1	1	1	1	1	1	1	1
D4Mit17	33.96		1	2	1	1	1	1	1	2	2	1
D4Mit27	42.13		3	3	2	2	3	1	1	1	1	2
D4Mit31	50.04		3	2	3	3	2	1	1	1	1	1
D4Mit308	57.66		4	1	2	2	3	3	3	3	3	4
D4Mit203	63.26		3	1	1	1	2	2	3	2	3	2
D4Mit251	69.05		2	1	1	1	2	2	2	2	2	2
D4Mit354	76.12		3	2	2	2	3	1	3	2	1	2
D4Mit42	82.64		1	3	2	2	1	1	1	3	2	3
D5Mit145	3.37		2	4	3	3	4	4	2	2	4	1
D5Mit73	12.88		2	2	1	1	2	2	2	2	2	2
D5Mit352	18.4		3	1	1	1	3	2	2	2	2	3
D5Mit233	28.55		3	1	1	1	2	3	2	3	3	3
D5Mit307	40.31		1	1	1	1	3	2	3	2	2	2
D5Mit277	51.70		2	1	1	1	1	2	2	1	2	1
D5Mit158	55.99		2	1	1	1	1	2	2	1	2	2
D5Mit161	65.34		2	1	3	3	3	2	2	2	2	4
D5Mit99	79.51		2	1	2	2	2	1	3	3	3	3
D5Mit43	84.68		2	1	2	2	1	1	3	1	1	2

Figure 1. Continued.
Agarose Gel images of PCR bands and their interpretation.

			<div><div><div>/</div><div>N</div><div>J</div></div><div><div>6</div><div>/</div><div>J</div></div><div><div>S</div><div>v</div><div>J</div></div><div><div>9</div><div>S</div><div>2</div></div><div><div>2</div><div>/</div><div>J</div></div><div><div>A</div><div>/</div><div>J</div></div><div><div>L</div><div>/</div><div>J</div></div><div><div>b</div><div>/</div><div>C</div></div><div><div>/</div><div>H</div><div>e</div></div><div><div>N</div><div>O</div><div>D</div></div></div>									
Marker	cM	Gel Image										
D6Mit138	1.81		3	1	3	2	2	2	3	3	4	3
D6Mit296	2.25		1	1	1	1	2	1	1	1	1	2
D6Mit116	11.5		2	2	2	2	1	2	1	1	2	2
D6Mit274	23.70		1	2	2	2	1	1	1	1	1	1
D6Mit8	35.94		3	1	2	2	1	1	3	3	2	3
D6Mit100	41.03		2	1	1	1	1	1	1	1	1	2
D6Mit36	48.93		1	3	2	2	3	2	1	2	2	3
D6Mit10	52.75		1	2	2	2	3	1	1	2	1	1
D6Mit333	60.55		2	2	1	1	2	2	1	2	2	2
D6Mit302	67.66		2	3	1	1	3	3	2	3	3	2
D6Mit14	77.64		3	2	2	2	1	2	1	2	1	3
D7Mit21	1.91		1	1	2	2	2	2	1	1	2	1
D7Mit152	2.71		1	2	3	3	1	1	1	1	1	1
D7Mit57	9.94		1	2	2	2	1	2	2	2	2	2
D7Mit267	17.09		3	2	3	3	1	2	1	1	3	2
D7Mit158	29.82		1	1	2	1	1	1	2	1	1	1
D7Mit317	41.50		1	4	2	2	3	2	2	3	3	2
D7Mit194	44.83		1	1	1	2	2	2	1	1	2	1
D7Mit220	55.69		2	3	2	2	2	2	1	2	2	2
D7Mit101	69.01		2	2	2	2	2	2	1	1	2	3
D7Mit223	88.85		1	1	2	2	1	1	1	1	1	1
D8Mit143	11.59		2	2	2	2	1	2	2	2	1	2
D8Mit63	19.02		2	1	2	2	2	3	3	3	3	3
D8Mit178	34.43		3	1	3	2	3	3	3	3	3	1
D8Mit45	42.16		3	1	2	1	3	3	1	2	3	2
D8Mit80	43.06		1	1	2	2	1	1	1	2	3	2
D8Mit211	52.0		2	1	2	2	2	2	2	2	2	2
D8Mit215	62.63		1	3	3	3	2	2	2	3	1	2
D8Mit56	76.14		1	1	3	3	3	3	1	1	1	2
D9Mit249	2.46		1	1	2	2	1	1	1	1	1	1
D9Mit250	2.46		1	1	2	2	1	1	1	1	1	1
D9Mit89	14.79		3	2	2	2	2	1	3	2	2	2
D9Mit285	21.42		1	3	2	2	4	1	2	1	1	1
D9Mit336	35.39		1	2	2	2	2	3	1	2	2	3
D9Mit8	42.65		3	2	1	1	3	3	3	3	3	3
D9Mit347	55.11		1	3	3	3	4	4	2	3	3	3
D9Mit18	71.49		2	1	1	1	2	2	1	2	2	1

Figure 1. Continued.
Agarose Gel images of PCR bands and their interpretation.

Marker	cM	Gel Image	F V B / N J	C 5 7 B L 6	1 2 9 X 1	1 2 9 S 2	D B A 2 / J	C C B A / J	S J L / J	B a l b / C	C 3 H / H e	N O D
D10Mit28	3.04		2	1	1	2	2	1	1	1	1	2
D10Mit213	9.75		1	2	1	1	1	2	1	1	1	3
D10Mit108	22.89		2	2	1	1	2	2	2	2	2	2
D10Mit20	34.83		2	1	2	1	2	1	1	2	1	2
D10Mit230	45.28		2	1	2	1	2	1	1	2	1	2
D10Mit96	51.19		1	2	2	2	2	1	2	2	1	1
D10Mit178	51.42		1	2	2	2	2	1	2	2	1	2
D10Mit233	61.58		2	1	2	1	2	1	2	1	2	1
D10Mit271	72.31		1	2	2	2	2	1	1	2	1	3
D11Mit1	2.51		1	1	2	1	1	1	1	2	1	1
D11Mit226	5.64		1	2	1	1	1	1	1	1	1	1
D11Mit151	15.29		3	2	3	3	2	1	3	3	2	3
D11Mit217	23.1		1	3	3	3	1	3	3	3	3	2
D11Mit349	33.29		1	2	1	1	1	1	1	2	1	1
D11Mit298	42.76		4	2	3	3	1	2	2	2	2	4
D11Mit39	52.22		1	2	1	1	1	1	1	1	1	1
D11Mit70	58.9		2	3	3	3	3	1	3	3	1	3
D11Mit54	59.82		1	2	1	1	1	1	1	1	1	1
D11Mit132	62.92		2	1	2	2	2	1	1	1	2	2
D11Mit333	71.83		2	3	2	2	2	1	3	2	2	2
D11Mit303	82.9		2	1	3	3	2	3	4	3	3	2
D12Mit12	8.49		2	1	2	2	2	3	2	4	3	4
D12Mit2	18.94		2	1	1	1	2	1	2	1	1	2
D12Mit114	28.94		3	1	2	2	3	3	1	2	3	1
D12Mit149	36.85		1	2	3	4	3	3	3	3	3	1
D12Mit117	43.32		1	2	2	2	2	2	2	2	2	1
D12Mit30	50.43		1	1	2	2	1	1	1	1	2	1
D12Mit99	52.9		1	1	2	2	2	2	1	2	1	1
D12nDs2	62.22		3	4	3	3	2	1	4	3	2	4
D12Mit263	62.11		3	1	2	2	3	3	1	2	3	1
D13Mit16	7.26		1	3	1	1	3	3	1	1	3	2
D13Mit137	18.87		1	1	2	2	1	1	1	1	1	1
D13Mit63	21		2	1	2	2	2	2	2	2	2	2
D13Mit186	31.87		3	1	1	1	1	1	1	2	1	1
D13Mit282	32.59		1	3	3	3	3	3	3	3	3	2
D13Mit191	45.05		2	1	3	3	2	2	2	2	2	2
D13Mit51	56.45		1	2	1	1	2	1	1	1	1	2
D13Mit213	59.69		2	1	2	2	2	2	2	2	2	1
D13Mit78	67.21		1	2	2	2	2	2	1	1	2	1

Figure 1. Continued.
Agarose Gel images of PCR
bands and their interpretation.

Figure 1. Continued. Agarose Gel images of PCR bands and their interpretation.			I									
			C	2								
Marker	cM	Gel Image	F	7	X		D			B	C	
			V	B	1	1	B	C	S	a	3	
			B	L	/	2	A	B	J	I	H	
			/	6	S	9	2	A	L	b	/	
			N	/	v	S	/	/	/	/	H	
			J	J	J	2	J	J	J	C	e	
											O	
											D	
D14Mit179	4.92		1	1	1	2	1	1	1	1	1	1
D14Mit126	11.94		1	2	3	2	2	2	2	3	3	2
D14Mit60	24.6		1	3	2	2	1	2	1	2	2	1
D14Mit259	27.65		3	1	2	2	2	2	3	2	3	2
D14Mit39	35.69		2	2	1	1	1	2	2	2	1	2
D14Mit68	37.61		2	1	1	1	1	1	2	1	1	1
D14Mit194	45.96		1	2	2	2	1	1	2	1	1	1
D14Mit95	57.2		1	1	2	2	2	1	1	1	1	2
D14Mit97	62.20		1	3	2	2	3	3	1	3	3	2
D15Mit13	1.84		3	2	1	1	1	3	3	3	3	4
D15Mit265	6.08		3	3	3	3	1	1	2	1	1	2
D15Mit138	15.68		2	2	4	4	3	3	3	3	1	3
D15Mit154	16.82		1	2	1	1	2	2	3	3	1	2
D15Mit128	25.58		2	1	1	1	1	2	1	1	2	1
D15Mit92	32.19		2	3	2	2	1	2	2	2	2	3
D15Mit80	38.02		1	2	2	2	2	1	3	1	1	1
D15Mit34	45.31		1	1	1	1	2	1	3	3	1	1
D15Mit244	48.65		2	1	1	2	2	2	2	1	2	2
D15Mit16	58.05		2	1	1	1	3	1	3	2	2	3
D16Mit131	3.41		2	1	3	3	3	3	3	3	3	3
D16Mit136	27.82		1	2	3	3	2	3	2	3	4	3
D16Mit30	33.01		1	2	2	2	2	3	2	2	1	2
D16Mit64	34.22		1	2	2	2	2	3	2	2	2	2
D16Mit140	40.30		2	1	2	2	2	2	2	2	2	2
D16Mit19	45.36		2	1	2	2	2	2	2	1	1	2
D16Mit70	48.81		1	2	1	1	1	1	1	1	1	1
D16Mit71	57.06		2	1	2	1	1	1	1	1	1	1
D17Mit164	2.11		2	3	2	2	3	3	3	3	3	1
D17Mit133	12.53		1	3	1	1	2	1	1	1	1	2
D17Mit51	19.74		1	2	2	2	2	1	2	2	1	2
D17Mit68	23.55		2	1	2	2	2	2	4	3	2	4
D17Mit20	29.73		1	2	2	2	1	1	2	1	1	1
D17Mit205	39.3		2	3	2	1	2	2	2	3	2	3
D17Mit1002	50.97		1	1	1	2	1	1	1	1	1	1
D17Mit123	60.67		3	1	2	4	4	3	4	3	4	2

Figure 1. Continued.
Agarose Gel images of PCR bands and their interpretation.

Marker	cM	Gel Image	<div> <div>1</div> <div>C 2</div> <div>5 9</div> <div>F 7 X D B C S B C</div> <div>V B 1 1 B C S a 3</div> <div>B L / 2 A B J l H</div> <div>/ 6 S 9 2 A L b / N</div> <div>N / v S / / / / H O</div> <div>J J J 2 J J J C e D</div> </div>									
D18Mit222	8.08		2	2	2	2	2	1	2	1	1	1
D18Mit226	18.18		2	1	2	1	1	1	2	2	2	2
D18Mit177	21.39		1	2	3	3	3	3	1	3	2	3
D18Mit91	29.67		3	2	3	3	1	1	4	3	2	2
D18Mit124	32.15		1	3	2	2	2	3	2	3	2	2
D18Mit184	39.70		3	4	1	1	1	1	3	2	1	3
D18Mit186	45.63		1	2	1	3	3	3	3	2	3	2
D18Mit7	51.92		2	1	4	4	3	3	1	1	3	1
D18Mit129	53.28		2	1	2	2	1	1	1	1	1	1
D19Mit68	3.38		1	2	2	2	2	2	1	1	1	1
D19Mit60	13.9		1	2	2	2	2	2	1	2	2	1
D19Mit45	16.14		3	2	2	2	1	1	3	2	1	1
D19Mit30	21.34		1	3	1	3	3	3	3	2	3	2
D19Mit88	32.23		1	2	2	2	2	2	2	2	2	2
D19Mit19	34.08		1	2	2	1	1	1	1	1	1	1
D19Mit11	36.26		3	1	1	4	2	4	4	3	4	3
D19Mit91	40.53		3	2	2	3	3	1	3	3	1	1
D19Mit1	50.32		2	1	1	2	2	2	2	2	2	2
D19Mit71	56.28		2	1	1	2	1	1	2	2	1	2
DXMit136	4.23		1	1	2	2	1	1	1	1	1	1
DXMit48	25.51		1	2	1	1	2	2	1	1	2	1
DXMit62	34.6		2	2	1	1	1	2	1	2	2	1
DXMit110	35.53		2	1	2	2	2	2	1	1	1	1
DXMit170	45.87		2	1	2	3	1	3	1	3	2	1
DXMit179	53.17		1	1	1	1	1	1	1	1	1	2
DXMit79	68.46		1	1	2	1	1	2	1	2	2	1
DXMit153	73.26		1	1	2	1	1	1	1	2	1	2

Table 2: The markers that yielded PCR bands but could not be included in the panel due to their not meeting the set criteria

Category 1: Markers that yielded same size bands in all the strains

Marker	cM	Primer F	Primer R	Lab code
D2Mit74	103.4	CCAAGCTTGCAGTTTGTAGC	AGGTGTTATTGAGCCCTGTATAGC	2-6A
D10Mit88	28.64	AAGATGAGAAGATAACATGTCAGGC	TTCTGAATTAAGTTCATCTGAACCC	10-2C
D11Mit226	5.64	AGGTGAACCTCTTTGAAGTTTGTG	AAAGGAGTGACTGAGAAAGACACC	11-1
D11Mit60	42.86	AGAGAGGCAAAAATTCGAAGC	CTTCCTGATGGTAGGATTTAGGC	11-60
D16Mit51	53.78	CCTCAGGTCACTCAGGATTTAA	CCTGTTCCACCTCTCCACAT	16-5
DXMit57	16.24	AGTAGCAAGTAGACTCTCAAAGAGGG	TCTGGCATACTGGGCACT	X-57
DXMit81	20.59	GAGGAGCATCAACCTTCTCG	GAGGTGGGGAGAAACAGAGG	X-81
DXMit63	41.51	TTATAAATTAGTGTTACCACATGCAGC	AACATTTTTTCTAGCATGTGTG	X-63
DXMit186	76.75	ATCAATGCATAGTATTTGGGCC	AATTTGCTACTGCGGGTAGG	X-186

Category 2: Markers that failed to amplify a PCR band in some strains

Marker	cM	Primer F	Primer R	Lab code
D1Mit375	23.18	TAAATCCATAGATGATAGATCAGTGTG	GTGGAAAAAAACCTAAGACACC	1-2C
D1Mit76	33.31	ACAAAGGAACTAAACAGACTCGG	CTCCCTCAAATACATCTTTGGC	1-3C
D1Mit308	50.78	GAGGCTATGAGTCAAATGGACC	TTTATGAGGTGCTGAGATGCA	1-5
D2Mit448	65.66	TACTGTTTGCATTGAGTGCG	AAAAAGTAATGGTTGGGGCC	2-4B
D3Mit278	32.59	TCTAATATTGGAAAATGAATTTCTCTC	TATGCCACATGCACACC	3-3P
D5Mit309	42.22	TAGAGCCTATTTCAAACCCCC	GTTGCATCCATAGCAAGCAA	5-4
D10Mit49	1.91	GGAATTTACTGGAATACAACCC	GTGGGCATTTGCACTGTG	10-1
D11Mit168	78.74	CAGGGATTTGACTTTTAACCTCC	GAAATGGCTCTACAACCTCC	11-168
D14Mit132	6.03	GAACAGCACCATCCACACAC	GTGGGGTTATATGCAGATACTCG	14-1
D15Mit147	46.85	GATGTGTGAAAAATTTGTTTCTTG	GTCTCAAAGGAATAAGAAAGAGATGG	15-5
D17Mit3	34.9	GATCTTTTCTTATTCTGGTT	GCAAAGTCATGTACTCTGAG	17-3B
D17Mit39	45.64	CCTCTGAGGAGTAACCAAGCC	CACAGAGTTCTACCTCCAACCC	17-5
D18Mit48	51.89	TTGCACTCACAGGGCACAT	TCAGAGTTTCCAGAAGACACCA	18-6
D18Mit25	59.05	CTGGAAATAAAACCTGGGCA	TTTAGCCTAACTGAGTTCAGACC	18-7
D19Mit31	8.84	CAATACAGAGTAATGATTGCCTGA	TTACATTTTGGGATGCTCA	19-1B
D19Mit41	13.18	AGCCCTCCACCCAGTTTC	TCTGGGGAAAAAGGATGAGA	19-2
DXMit124a	4.25	AAGGAGAAGTAGAAGAAAGAAAGAGG	GGTGTAGCCTCAAAAAAGATGG	X-124a
DXMit68	29.49	TCCTTTGGCCTCTGCATAT	TGTTCTTACAATGAGCCTCATAGG	X-68
DXMit114	42.82	ATGGCATCCACAGTACCACA	GTAATCAATTTGTGAATAAGGAAGC	X-114
DXMit95	45.86	CTGTCAATCTAATCTCTATGTCTGTG	CTTCTGGGTGGCAGTG	X-95
DXMit155	67.49	TGTGCAACCACCATTTC	AGGATCTGAGTGCCCAACC	X-155

Category 3: Markers that amplified multiple PCR products in most (or all) strains

Marker	cM	Primer F	Primer R	Lab code
D11Mit86	32.13	TTGACATTGTGACAAAGACTTTCA	AAGGCATCATGAGGTTTTTAGTG	11-3
D12Mit118	44.93	CATCTTCAATAAAATGGAGATGTACA	CGCTTCCCTTCATGTACTAGC	12-4
D13Mit107	50.2	CAACTGAGCCACATTCTAGC	CAGCCACAGGATGATGAGG	13-4
DXMit98	72.38	GAGAGCAGGACTATGACTTG	ACACACCAGTTCGCACACAT	X-98
DXMit222	78.31	TTGGTTTGGGGTTTTTTTTG	ATTCTGATAATGTCTTCTGGACA	X-222

Table S1: Markers that did not amplify reliable PCR bands in most/all samples

Marker	cM	Primer F	Primer R
D1Mit22	30.44	TCTGTTCCCTCTACACACATGC	CTACCATGCTTACCTAGGTCCTG
D1Mit136	49.26	TAGCCCTACACACTGTAGAAATGC	TGAACACAAAGTAGTAAATGCGTG
D1Mit87	50.93	CACATGCAGTGGCACATATT	TTTTACCATGGGGTTATTCAAGG
D1Mit285	51.56	GCCCCATATGTATGCACTGCT	GCACCCAAGGTTGTCCTCTA
D1Mit309	52.14	GTAGAGAATGATAAAGGACATCCTCC	GGGAGGCTGAGATCTGTCAG
D1Mit94	54.7	CGACTTCCCTTGATGTCCAT	TTTGTGTTGTGCAGTCTGTCTG
D1Mit288	63.51	TGAAATTGGGAACCTCACGTG	CCATTGTTTTGCTATCCTTTCC
D1Mit151	84.49	TTGCCCATCCTAGTGAGAGG	CACCATACGTTTACAGATAGACAGA
D1Mit37	84.93	ACAGGACTTCTTACTCAAACCACC	TTCTTTTGGCCTCTTTGGG
D2Mit354	3.11	CGCTGACAACAATCTCCATG	CCCACAGTGTCTTTCTACAAGC
D2Mit267	12.34	TGTCTAAGGGATAGGCCACG	AGCACATATGCTCGCACTTG
D2Mit464	12.9	GAATAGGGGAGTTAGAGGGGG	GAAACTCTTTTTTTTCTTCAGTGTG
D2Mit324	33.33	TTTCTATTCAACCCTAGGGCT	TTCTTCTGTCTCTGTTTCTGTCTC
D2Mit61	34.9	AAAGTCAACTGCTTTCAGTTACCC	CACAGAAGTGCCCTTGCTA
D2Mit141	51.62	GCACTCTGCCAACTGAATCA	TGTGTGTGCACAAGCATGAG
D2Mit101	57.97	ATAATTCCTGATTTGCTGTTTGTG	ACATGAAGCCTAGAGGGTGC
D2Mit422	70.07	TGAGTGCAAAGTGCCTCAAC	AAGTTGTCAAAGTGTTAGACAAGGG
D2Mit59	75.99	TCCAGAGTGCTTATATGA	AGAAAANACAGAAAAATCAC
D2Mit227	85.27	TCTGTCTGCCTGCCTATGTG	TGAGATCTTGTCTCAAGAAACAGC
D2Mit265	97.97	AATAATAATCAAGGTTGTCATTGAACC	TAGTCAAAATTCTTTGTGTGTTGC
D3Mit178	14.27	GTGTACATGCACATGAGTTTGC	CAGAGAGTGCTAGGCTCAAACA
D3Mit63	19.65	CTATGGACTTGTGACATAGGAGTG	ACCACAAGATGGAGTAAGAGTTCA
D3Mit224	20.54	ATTGCTATACCAAGTAGACAAACACA	ACAAGGAGGCCGCACTAAC
D3Mit10	43.03	CTGGCTTGGTGGAAGTCCT	CCTAAGCCAGCTACCACCAC
D3Mit215 F	50.18	TAAACATCTAGAAGATGCTGCAGG	CTGCATGGCCAGGACTAGTT
D3Mit42	52.94	TGACCTCCAGAGAGTCTTCCA	TAAGCAGCTGAGACTCAAGTGG
D3Mit13	54.68	TTTCTGCATTATGTGGGCTT	AACCACAGATGACAATTGAA
D3Mit17	67.53	CATGGCTCCATGGTTCTTG	CCACGGAGAACAACCTGAAGA
D3Mit18	72.36	GAACAGTTCCTCAGGTCCTCA	CTGCCTTTAAATTCTGTCACCC
D3Mit147	74.02	TCTGCCTCTGTTAGATAGATATCCG	TTGTTTCATCTATCCTCTGAAGTTCC
D3Mit201	80.58	ACCATCATTAATAAAGGTTTGACA	GTTTACCTTTGGACTCTACACATGC
D4Mit227	4.43	CTCAGACATGATTTTTTCCAAGG	GCAGTTAAACTGTACTTTCTGTAAACA
D4Mit214	24.16	CCCTACAAGTTGCTCTCTGACC	CACCAAGTGTGAGAGCAGGA
D4Mit270	24.49	GGGACCTGAGGTGGGGAG	GCCCTGTTTTAAACTGACCC
D4Mit38	24.82	AGAGACTGGTGGGAATCCT	TGACTAGACAGTCTTCCACTGAGC
D4Mit27	42.13	GCACGGTAGTTTTTCCAGGA	TGGTGGGCAGGCAATAGT
D5Mit73	12.88	GTTTGAGAGGTCCTGAAAGCA	TTTCCATTTACATACATTTGTGCA
D5Mit229	17.27	TGGTGAGGTCAAACCATAGC	AATTGACCTCTGCCTTTATATGC
D5Mit201	39.55	GAGGACTCCTTCGATTTCCC	TTCCTAAGCAGGAACCTGACCA
D5Mit113	41.59	ACAGTATTTCTTTTTCCAAGTGTG	CAAAGACTCTAGGTGTGACCCC
D5Mit219	74.84	TGCCTTGTTTCTATCCACC	GACTATATAGTGAGATCCGGGCC
D5Mit99	79.51	CAGAAAAGAGAAAACGGAGGG	TTCTGCTGCCTGAAGTTTT
D5Mit144	89.64	GAATGGCCCCATAGGTTCTT	CCAGTGACACACTTCTCCA

Table S1 Continued: Markers that did not amplify reliable PCR bands in most/all samples

Marker	cM	Primer F	Primer R
D6Mit116	11.5	ACATTTCTTTGTGAGGTTCTTG	CAGGTTTTTTGAAAGACACTCTTG
D6Mit180	12.36	TGCTTTGAAGATGGCAGATG	ACCCACGCAATTCCTAAG
D6Mit48.3	12.75	CACTCTGTGTAAGGAGCTTTTCTG	GCGTCATATGATAATTCGCG
D6Mit117	19.76	CTAGCAACCATATGCCTGCA	AAAAGTATAATTACTGGTAGCCCCG
D6Mit123	27.76	GGAAGGAGCAGGTCCAATAC	CTCCCAACCACCAAGACCTA
D6Mit300	58.07	TGAGATCCTACCTCATTACACTTCC	TTGGTTCATAGCCTTTCATGG
D6Mit12 F	59.14	CCACATCCATGTAAAAGCTG	TGGTTCAATGAAAGTTGCCA
D6Mit24	59.17	CTGTGGAGCTAGTTCAGGC	TTTGAAACCAGAATGGCTCC
D7Mit21	1.91	GGGTTGAACCTTACAGGGGT	ATCAAACCAGCCCAAGTGAC
D7Mit243 F	3.16	TTCCCTTCACTCATGATGTAAGA	GGCTGCCATACTGTTGAACA
D7Mit57	9.94	TTCCCTCTAGAACTCTGACCTCC	AGTTCAGAGCCGAGACTAGGC
D7Mit341	10.32	GGTTCACCAGCACTGCACTA	GGCAGATGTAGCAGTGCAAA
D7Mit80	26.29	AGCAACACAGATTGGACTTGG	CATCCTTCTGACTCTTCCCC
D7Mit69	31.19	CCCACCAGAGATCACCAGT	CACAAATGAAGGCTGAAAGCA
D7Mit89	34.6	CCTCTGTTACAGTCAGTGGTTAAA	GGCTTTTATTTCCATATATGGGG
D7Mit350	47.43	TCTGCATCTCACTGTCCCAG	ATCTACAAATGAGTTTCTAAGGACTGC
D7Mit43	73.19	GATATAGGGTGTCTTCTCAATCA	AGACCTCTCTACCAGAAGCC
D7Mit359	76.03	ATACAACTGAGCAGGTTGTAGTTATG	AAATACATGGCTTTTTCTTTAATCA
D7Mit108 F	76.89	GTGTTCTTTACCACAGGGGC	TCTCAGAGAATATGCTTCTTCCC
D7Mit207	77.13	GCTTATTACAAACAAGCCTGCC	GTCTTCTACACATGTGGAAACC
D7Mit304	77.48	GCCCTTCATATGCATGCC	GTGCATTTTAGTTTGTATTGCTGC
D7Mit332	77.87	TACCATCCCTAACTGGTTCTCTT	CTGCACACTCACATACATACTCATG
D7Mit166	78.36	TATATAAATCTATCAAGACACACCCCC	AGTATCTGATTTTGATCTCTGGCC
D7Mit287	78.67	TGAAGTATGACTGAAGAGGCTCC	AGGCTGATATTGGTGTAGTAGGG
D7Mit360	81.36	GTCCAGGGAATCAATCCAAT	CTGTGTGTGTCTCTGTATCTGTGTG
D7Mit290	81.59	AAGAAATGGCATCTCTCTCAGC	AAGTGTAAGAATTTAATGCATACACACA
D7Mit12	81.85	GCTGGGTTTATTCATTGCAA	TCCAGCTCATGGGTAGAAGA
D7Mit361	82.32	TACTACCAATGCTAAGCGCTACC	ACCAGAGTTTGTGCCCATTC
D7Mit242	82.95	ACACAGGCCCAACTCGTAAG	CAGAATGGACACCCAGC
D8Mit155	2.14	TTGGACAGGGAAAATTCTGC	TGAGGACTTGCTTTAAGAGTACTCC
D8Mit141	5.73	TCCAATATTGACCTCTGGCC	CTGTTTATAAACTCACCCAGGTTCT
D8Mit143	11.59	AGCCTGAGGTTATGTTTTTTCG	GGCCCTCAGGTTTTCTCTCT
D8Mit356	14.77	AAAATGTCTCTACATTCTACATATGCA	GTCCAGTCCCAGTACCCAGT
D8Mit4	18.89	CCAATCATCCCCAAAGGTA	GTATGTTCAAGGCTGGGCAT
D8Mit24	21.16	ACCTCCACACGATTACAGTGG	ATGATGTCACCAGCAACTTCC
D8Mit190	22.16	CTTTGTTGCTGTTTCATTCTGG	AGTCATATACAAGGTCAACCTGAGC
D8Mit225	23.05	CCCTCTTCTTCCCTTCCACT	TTTGTTGTTGCTTGCTTTGG
D8Mit97	23.05	CTGTGCAAAGTTGCTAAAACAC	ATCTTTTAACTTAGGTGGAGAAAAACC
D8Mit5	24.7	TCCCTTTTCCCTGTGCTATG	GCCGTTCAATTAACCTTCA
D8Mit229	25.03	GCCTTCACTTAAGGCATAGGC	AAGAGCAAACCTGCCAGTAAGG
D8Mit175	25.86	GCAAGAAAAGTCTCAGATGC	GAGGGGTGCTTTGGCTC
D8Mit6	27.69	CAGGCAGCTTGCTAGGACTT	TACTGCCTTTAGCCCAGTGG
D8Mit205	28.85	TTGACTGCATTGACATTCTCC	GAGAAGAAGGCACAGAGTTTGG
D8Mit67	28.91	TATTGTGAGTGGAGGAATATCTGC	GATTTGTGAGTATGAGTGCATGC

Table S1 Continued: Markers that did not amplify reliable PCR bands in most/all samples

Marker	cM	Primer F	Primer R
D8Mit155	2.14	TTGGACAGGGAAAATTCTGC	TGAGGACTTGCTTTAAGAGTACTCC
D8Mit230	29.06	CAGGTTCCAAAAGTCTTTCACC	TTAGAAAGGTACTGATTGCTCTGC
D8Mit248	44.99	ATCCCTCAAGCAGTACCCCT	AGCAGAGGACCACACCTTACA
D8Mit201	69.04	AATTGAGGAAGACATCTGACATTG	GCCAGATCTCCTTACCAGACC
D8Mit55	70.08	GGGCACAGTGATGTTACACAC	TTGGGAATCTGGTTACTTCTCG
D8Mit323	72.73	TGTTAGTGGAGGGGTGGTGT	GCAAGATTTGAAACCAGGGA
D8Mit92	74.62	CTCAGGCTATCTTGGACATGC	TGGCTCACATCTGTGCTTTC
D8Mit156	76.09	ATTCTGCAATGCAAATACATACA	AACCAATGGCCTCCAGAAG
D8Mit156	76.09	ATTCTGCAATGCAAATACATACA	AACCAATGGCCTCCAGAAG
D9Mit64	14.35	TTCACCAAACCTTATCTTACTCCA	TGGAAGAAACAGTGTTGGGT
D9Mit229	26.83	CCCAGTAGCATGGAGAGAGG	TGTTTCCTTTGCTCTCAGCA
D9Mit36	52.64	GACCACAACCCATCTTGCTT	TCCAGTTTCAACTCTGGGCT
D9Mit350	61.01	GATCAGAGGTGCACGTTGTG	CTCCTGGTCCTTGTTTGTGTTG
D9Mit185	61.21	AGCACACATACATATACGTCCACC	TCGGCTCCTTTGCTTTTG
D9Mit81	65.15	TGGTGTGTGGTTTTTCAAATG	CCTGTGGAAGAAGTGGGAAA
D10Mit3	16.53	GTTGATAGTCCCACCTCACTCA	TGAGAAATTCCATCTGTGGC
D10Mit192	18.85	CCCAGCCATTGAACTGTTTT	CCAGTGTTCTAGGTTCCCGA
D10Mit156	19.29	ACCACAGATGAAAGAGAAATTGTG	TCACTATAATTTCCAATTGCTTTATTT
D10Mit148 F	23.24	TTCCAGACATACCATCACTATTCTG	TTAAGTAGACTCTGTTGCAGAGGTG
D10Mit58	25.54	TGTGTAATGTGCTTGTCTGTGG	ACAGAGCCTGCAGGTGTTTT
D10Mit186	38.56	CCTCCATTTAAAAATAGTGTTGTGG	CCCTGACCTCCACATATGCT
D10Mit10	46.85	CCAGTCTCAAAACAACAACAAC	TTGCACCTAGATTGCCTGA
D11Mit72	2.94	TGAGTGGCATCCCTAATTCC	GTGTGTGTCACTATTCTTGTTCC
D11Mit229	15.63	TGTTTGCTTGGTTTTGTTTTG	ATCCCGTTTTTACATCAGACA
D11Mit217	23.1	ACTGGAAAATATGTTTTAAACCTCTG	AAATGGGATTCTGCAAAAACC
D11Mit20	27.23	CCTGTCCAGGTTTGAGAGGA	CTTGGGAGCCTCTTCGGT
D11Mit207	32.13	TCTTTAAAAAGATTCATTTTGTAGGG	TCCCTCCTCTGGCTTCTACA
D11Mit242	39.47	GAAGCCAGCAAGAAAAATGC	CTGTCTGGTAGTGCAGCCAA
D11Mit5	40.59	TTCTGTGAGCCTGGAGGAGT	TACAGGACTAGTTTCCATTGCGG
D11Mit4	41.87	CAGTGGGTCACTAGTACAGCA	AAGCCAGCCCAGTCTTCATA
D11Mit278	42.51	AGGTCTCTCCAGATCCTCTC	TATCAACCACTCCAGGACAGG
D11Mit7	45.92	AGGGTATTCTCTAGCCTCCACAC	TTTGAGGCAAGATGTCATGTATG
D11Mit34	46.74	CGCCAGCCAGGAATACTTAG	AGTAAAATGAGCCCTTTCCTCC
D11Mit54	59.82	AGGCTGGTGGCTAGTGTC	AAGTCTTGCGCTGCATCTTT
D11Mit12	78.74	AGGGTTATGCTCTTGGCTGC	GATTTTCCTAGGCTGGCTGG
D11Mit214	80.6	CATACAGCCTTCAACAATGACA	ACTGCATACATGTGCACTCATG
D11Mit153	29.93	ACATAACCCAGCCCTCCC	TCTGTGTTTTTCTCCAAAGGTG
D12Mit165	45.97	ACACTGGATCTTTTGCCTAGG	CATAGATCACATTAAGCAACAGCA
D12Mit120	48.78	AACATACCTTCACTGATGAGTTTCC	TTGATAACATCGAGGCTAGTATATGG
D12Mit30	50.43	TATGTGACTGCAATCCCAGC	ATGAACACATCATGCCCAGA
D12Mit231	50.87	GAGTGGATAATGAAAATGTGGTG	CCTGACATTTTTATGATTTTATTTTC
D12Mit102	53.88	TGCAGAAAGACAGACAGAAAGTG	TCTGCCTCATTGGATCACTG

Table S1 Continued: Markers that did not amplify reliable PCR bands in most/all samples

Marker	cM	Primer F	Primer R
D13Mit235	2.36	ATGTAGAGAAAGGCCACTGATTG	GGGGGGATAGCATTTGAAAT
D13Mit152	3.25	TGGCTGTGTCTTGCTTGTTT	AAGATATTCTGTTATAGCCATAATGGG
D13Mit60	14.4	CTCACAGGAATGGGGCTCT	TGTGTAGACAAGCAGAATCACTTC
D13Mit275	14.5	TTAGCAAGGGAACAGAGAGAGG	CAATCAAGGTATCCCTGTCTCC
D13Mit38	17.62	TGACATTTGCCTGTGGCAT	CACTCAATATTCTCCAGGGAGG
D13Mit88	18.48	ACTGATGGCTCATGAGACCC	AAAATTAATAGGAACTGCAAGGG
D13Mit88	18.48	ACTGATGGCTCATGAGACCC	AAAATTAATAGGAACTGCAAGGG
D13Mit61	20.16	TGCTCCAATACAACAAGGTCC	CCAGCCAAGGTGTGTTGAC
D13Mit62	20.43	CTATGGAGTTTTTGGTCATTTGG	TCAAAGAACCCTTTCCTCACC
D13Mit19	21.6	GGTGAGTTGTGTAATGATGGACA	AGCAACAGGGGCTACTAAACACA
D13Mit19	21.6	GGTGAGTTGTGTAATGATGGACA	AGCAACAGGGGCTACTAAACACA
D13Mit155	23.04	CAAGCCTGCTCACATTCTGA	TCATTGTTTTCCCCTGAAGC
D13Mit211	56.16	TTTGCAAAAAGACAAGACTTATTGC	AAAAGACATCCAGTTCTCAATGG
D13Mit51	56.45	TCCTGCAAAAAGTGGAGCC	TGGAACAAGCTCTTGGAGG
D14Mit253	14.55	GATGGGTGAACATCTGAGCA	GAGCATCTTCTTTCATACATGGG
D14Mit209	15.06	ATGATGGATGGACAGACAGATG	GTGGGTGGATGCTTCCAG
D14Mit113	31.51	TGCACAGGTTTTCCAATTTG	TGCTGTCTCTCCCAAGC
D14Mit65	31.62	TCACATGCCACAACCTTAGCC	CAACCCACACACACCTCAGT
D14Mit178	65.26	GCATGTACCAAGTGTGTGGGT	GCAGTACTAACTGGTTCACTGGG
D14Mit14	16.80	GCACATTCCAAAACACATGC	GGGATGGTGTCAATCAATCC
D15Mit174	1.84	GTGTTCTTGTTGTTTCTTTTTGG	GGTCAAGGAGACACTCAGTTCC
D15Mit13	1.84	GGAGACAAAATGAACTCCTGG	TTGTAAGACAAGCATAGCTCAACA
D15Mit80	3.82	TGAAGTCATCTTTCAATTTTCTCC	CGAAGATGCCTGCCAATAT
D15Mit18	11.93	GGGCTAATTGATAAATGATTAGTGC	CCCAATTCAGGGTTTCTAACC
D15Mit168	26.13	AAGGCAGGTAAACTGCAAGG	AGTGAATAAGAGACTTGGAGGCC
D15Mit39	52.36	CATCAAGGGGATAGAGAGTTGG	TGGGAATTGAGCCTGAGTCT
D16Mit1000	4.19	GGTAAAGCCCAAGTGGGG	AGAGAGATGGGCATGCACTT
D16Mit143	11.23	GAATCCCTCCTCTGGACCTC	ACAGTTTTTGTAGGCTACGTATCG
D16Mit1	12.89	CGCCCTCTAAGGTGACTCAG	AGAGAGGGGTTATGGGGTTG
D16Mit98	13.5	TCTCACGTCCTTCCCAAATC	ACCACCTGTAATGACTGCACC
D16Mit2	13.8	CCAATGCCCTCTTATGACCT	TTCTAGTGCCTCTACCCAG
D16Mit101	14.69	TTATGAAATGTTTTATCTTTTGGGG	CTCCAGATGTAGAAATTAAAATCTTGG
D16Mit146	14.69	ACCAATCTGGAGATATGTTACAAGG	TGGAAGACACATACTCTCTCTCA
D16Mit156	16.88	CAAAAAGCACTAAATGGACTTGG	AGACTCTTCTTCTTATTATTGTCACA
D16Mit133	18.04	TTGAATCCCTAGCACTCATGG	CCCATATTCTCTTGTCTCACTCA
D16Mit3	19.88	TCTAACGCCCTCTCTCTACC	CCAAATGTGATTGCACAAGG
D16Mit57	20.26	AAAAAATTTTAAACCATGTGAATGT	TGAAGTTTATTATGAGTTGAATCATGC
D16Mit36	21.34	CAGAGAATGATGAAAAGGACACC	ATGGACGTCTCAGCCTCTGT
D16Mit75	25.45	GAAAGTTCCAGGCTAGTTGGG	TAGGCCCTTCATTGCAAATC
D16Mit136	27.82	AGATAATCCCGTGAGAATAAAACC	TTGAGAAGTTTGGCCCTATAATGG
D16Mit71	57.06	TAGAAAATCTTCAAATAGGATCTGTTC	GAGCATTTCCCTTTTACCTGG
D17Mit113	8.14	TCTGTCTCCTCCGTAAGTGGG	GTCAATAAGTTCAATCACTGAACACA

Table S1 Continued: Markers that did not amplify reliable PCR bands in most/all samples

Marker	cM	Primer F	Primer R
D17Mit177	24.51	CCCTCTCATGATTTAGGACCC	AAAGGGCAGGTCACTCAAGA
D17Mit193	38.8	CACCTGTCCCCAAACAAAAC	GAAGATCTGGACATGCTCAGTG
D17Mit128	50.81	TGTCTGCCCCCTCTCATTC	TTAACTGTGTCTTGTATCCCATGG
D17Mit155	55.02	TGAGAAGGTTGGGTTTATATATTTAGG	CGATCATTTCTTGCAACCT
D18Mit66	2.42	AAACAAAACCAGAACAAATAGAGCA	TAAGTTCCTCTCATTTTCTGACCC
D18Mit31	5.8	AGCAGGAAGACCATGAGTATAAGG	TGCTGGCATGAGACAAGTTC
D18Mit68	11.96	GCGTGAGGGTTTTGTTTGT	AATACTTCCAGAACCTTAGACCCC
D18Mit22	13.41	TGATGGGATGTTTCTTGGGT	CACTGGATGACACAGCCTGT
D18Mit139	35.11	GTATTTGTATACAAGCAGGTATAGGGG	ATGTGCATGGACACATAAGCA
D18Mit4	57.53	ACTGTTGCTGGGGAATGG	CCAAGTTCAAAGCTGCTGG
D19Mit46	29.06	ACCCTGCCCTCTCTCTCC	GCACTCGCCAACTCAGGTAT
D19Mit12	30.30	TGCATATGTTTTCTCAAATGTGG	GTTCCCAACACACCTTTATACC
D19Mit25	45.05	ACATTTGAATATTGGAGAGGAAGC	TCTTGACCTAAACCTTAGTGGTG
DXMit89	4.73	TTCAGTTTTCCCTTCCCATG	AGGGAGTTACTGGGAGGGG
DXMit103	7.25	GTGGGAGTTGTCTTCACTGC	CTGACATCCTCCCATGGTCT
DXMit104	23.19	TAAACCTCCACACTCCTGGG	GCACAGGCCTGCTCCTATAC
DXMit83	27.99	TTCTAAACTGCTTTGTAAAGAAGGC	GTGAATGTAGGGAATGTCAAAGG
DXMit119	37.24	CTTTAACCATAATAATGGCCTTGC	GGGTTCTGTGATCGCAAGTT
DXMit10	68.46	GAATTACAGGCATGCGTCCT	TGTTTGACTGAGAGGATGCG
DXMit99	73.73	GGAATTCAACAAAAGGTATGTGC	ACCTGTCTCCTTCTCTCTCC
DXMit21	75.66	TGGGACAGAACCAACATCAA	TCAGTTAGGCGCTGAAGGTT